TEXAS ASM UNIVERSITY



No. 296

May, 1983

Subramanyam, V., Linder, K., and Srinivasan, P. R. 99Tc and 31P NMR Spectroscopy of a	Block, R. E. In vivo ³¹ P NMR Spectrum of a Rat Tumor	0
Cationic Technetium Complex for Myocardial Imaging	Springer, Jr., C. S., a Balschi, J., Yarmush, D. M., and Chu, S. C. Facts About Aqueous Shift Reagents for High Resolution Cationic NMR	1
Alkynyllithium, Rotational Barriers in Acyl Dihydropyridines	Cohen, J. S. NMR Studies of Perfused Cells; Mag. Res. Biol., Vol. 2	
Post-doctoral Position 4	Henriksson, U. and Wärnheim, T.	-
Williams, E. A. and Donahue, P. E. Carbon-13 NMR of PPO Copolymers 5	Simultaneous Measurements of Self-Diffusion Coefficients and Spin-Lattice Relaxation Times	5
Trska, P. and Hajek, M. Configuration and Conformation of 2-Alkyl-3,3,4-Trifluoroxetanes	Mulder, P. P. J., Pardoen, J. A., Erkelens, C., and Lugtenburg, J. ¹ H- and ¹³ C-NMR Spectra of 11- ¹³ C	
Schmitt, P., Wesener, J. R., and Günther, H., TANDEM-SEFT	All-trans-Retinal	7
Frye, J. S. Note of Caution in Using Br-79 to Set the Magic Angle for C-13 MAS Experiments	Murray, M. and Higgins, K. Measurement of H-H Couplings from ¹³ C Spectra of the AA'X Type	9
Kentgens, A. P. M., Scholle, K. F. M. G. J., and Veeman, W. S.	Smith, W. B. LIS on a Bidentate Ligand	1
²⁷ Al Magic Angle Spinning NMR of Zeolites13	Joseph-Nathan, P., Santillan, R. L., and Espiñeira, J.	
Wilde, J. A., Bleich, H. E., and Schwemm, R. E. 2D Magnetization Transfer in the Rotating Frame	Fluorine-19 NMR Measurements of Steroidal Drugs	2
Harris, R. and O'Connor, M. 2D Silicon-29 Spectra	Shapiro, B. L. 1983-84 Newsletter Subscriptions and Related Matters	15
Vincendon, M.		
NMR Investigation of the Cellulose Conformation in an Aprotic Solvent Solution		
Mutaw, M. and Kleps, R. A. UIMC + UICC = UIC; Getting Rid of Bad Vibes		
Keeler, J., Shaka, A. J., and Freeman, R. Practical Tips for Broadband Decoupling25		

A monthly collection of informal private letters from Laboratories of NMR. Information contained herein is solely for the use of the reader. Quotation is *not* permitted, except by direct arrangement with the author of the letter, and the material quoted *must* be referred to as a "Private Communication". Reference to the TAMU NMR Newsletter by name in the open literature is strictly forbidden.

These restrictions apply equally to both the actual Newsletter participant-recipients and to all others who are allowed open access to the Newsletter issues. Strict adherence to this policy is considered essential to the successful continuation of the Newsletter as an informal medium of exchange of NMR information.

WILMAD REFERENCE STANDARDS We surply Reference Standards to major instrument manufacturers and, in practically every case, our standards meet or exceed the instrument man facturer's specifications.

18

For Reference Standards not shown in the listings, please inquire di-

When ordering Wilmad Reference Star the tube size and don't forget to include t

Wilmad Reference Standards are sea high-resolution, Wilmad sample tubes ac otherwise indicated in the listing).

32.00

Res./Sens.

-				_			TUBE	STYLE	
	MAI	RInv	esti	da	ors			0	
	5mm	7"	2.57	3	A				
8mm	, 10mm, 12mm	7"	3.0						
	15mm	7"	4.0		4.				
	2.5mm	178mm	50mm		5mm O.D.	10-14-14 M-10-16-14-1 1-16-14-14-14-14-14-14-14-14-14-14-14-14-14-	98)	
- It	A All		LID.	TO			MINIST		-
ne		on a city to	To Week	tic			78		-
	A 10	,		10					
			1	WILMAI	REFERE	ENCE	STANDARDS		
		INC	115		AAA		PARAMETER		-
Page No.		IIIU	PHOT	8 /			SET-UP	5mm	
38	WGH-01 /	Acetaldehyde + 1	% TMS				Resolution	\$ 32.00	
		etald	TRE	MO	r ot	AA		32.00	
		6 A P IV		enz	-d, \$	W <u>*</u> 1 N	ерип	33.00	
	WGH-22	20% Acetaldehyd					Resolution	33.00 33.00	
		20% Acetaldehyd	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1		4		Resolution	33.00	
		Orthodichloroben:	A STATE OF THE PARTY OF THE PAR	ii Aceione-i	46		Resolution	32.00	
		30% Orthodichlor		% TMS in A	cetone-d.		Resolution	33.00	
		20% Orthodichlor		% TMS in /		W	Resolution	33.00	
	WGH-61 9	95% Orthodichlor	and the second second	% TMS			Resolution		
	WIRLTE	OR CALL	· LUR BI		IN BD-1	12	DA Vijon	33.00	
	WGH-69	0% Orthodichlor	obenzene n A	cetone-d,			Resolution	33.00	
	WE'H I	SEND+ III%	FinoVifile	SERV TO	FTHRN	MAIL	Resolution	31.00	
	WGH-64	0% Chlorotorm i	n Acetone-d,				Resolution	33.00	

Reference standards for use in NMR spectroscopy are an important and integral part of this investigative procedure. Our reference standards are sealed and degassed in ultra-precision, high resolution Wilmad sample tubes. We supply them to many of the major instrument manufacturers and, in practically every case, our standards meet or exceed the manufacturer's spec fications.

Our last published listing of our reference standards was more than four years old and, since that time, we have deleted and added some new standards based on usage requirements. In our Bulletin BD-102, we have regrouped therh and provided you with descriptions . . . and prices . . This bulletin is all you need to place your order for more than 130 items, including our new 2.5mm tubes.



WGH-36

Methyl Formate (Neat)

WILMAD GLASS COMPANY, INC.

Rt. 40 & Oak Road, Buena, New Jersey 08310, U.S.A. Phone: (609) 697-3000 • TWX 510-687-8911

TAMU NMR NEWSLETTER - ADVERTISERS

Bruker Instruments, Inc				see p. 10.
JEOL Analytical Instrument	s, Inc.			see p. (i) and outside back cover.
Nicolet Magnetics Corp				see inside back cover.
Sundance Supply Company .				see p. 23.
Varian Instrument Division				see p. 28.
Wilmad Glass Company, Inc.				see inside front cover.

TAMU NMR NEWSLETTER - SPONSORS

Abbot Laboratories
The British Petroleum Co., Ltd. (England)
Bruker Instruments, Inc.
Eastman Kodak Company
IBM Instruments, Inc.
JEOL U.S.A., Inc., Analytical Instruments Division
Dr. R. Kosfeld, FB 5 Physikalische Chemie, University of Duisburg,
D-4100 Duisburg 1, Germany
The Lilly Research Laboratories, Eli Lilly & Company
The Monsanto Company
Nicolet Magnetics Corp.
The Procter & Gamble Company, Miami Valley Labs
Programmed Test Sources, Inc.
Shell Development Company
Unilever Research
Union Carbide Corporation
Varian, Analytical Instrument Division

TAMU NMR NEWSLETTER - CONTRIBUTORS

Chemagnetics, Inc. E. I. du Pont de Nemours & Company Intermagnetics General Corporation The NMR Discussion Group of the U.K.


DEADLINE DATES: No. 297 6 June 1983 No. 298 4 July 1983

All Newsletter Correspondence, Etc., Should be Addressed To:

Dr. Bernard L. Shapiro Department of Chemistry Texas A&M University College Station, Texas 77843 U.S.A.

AUTHOR IN	DEX		TAMU NMR NEWSLETTER NO. 296, MAY 1983
a Balschi, J Bleich, H. E Block, R. E Chu, S. C Cohen, J. S Donahue, P. E		17 30 31 33	Mulder, P. P. J
Erkelens, C		37	Schmitt, P 9
Espiñeira, J Fraenkel, G			Scholle, K. F. M. G. J
Fraenkel, G		4	Shaka, A. J 25
Freeman, R			Shapiro, B. L
Frye, J. S Günther, H			Smith, W. B
Hajek, M		7	Srinivasan, P. R 1
Harris, R			Subramanyam, V 1 Trska, P
Higgins, K			Veeman, W. S
Joseph-Nathan, P.			Vincendon, M 21
Keeler, J Kentgens, A. P. M.			Wärnheim, T
Kleps, R. A		24	Wilde, J. A 17
Linder, K Lugtenburg, J			Williams, E. A 5 Yarmush, D. M 31

Variable Temperature CP-MAS with the GX Series FT NMR Spectrometers



 ^{13}C (50.1 MHz) VT/MAS spectra of hexamethylbenzene. a) and c) $^1\text{H-}^{13}\text{C}$ cross polarization. b) Bloch decay. The peak at \sim 90ppm is due to the Delrin rotor.



235 Birchwood Ave., Cranford, NJ 07016 (201) 272-8820

New England Nuclear®

Prof. B. L. Shapiro
Dept. of Chemistry
Texas A & M University
College Station
TX 77843

March 28, 1983

Sub: 45.05 MHz ⁹⁹Tc and 81.016 MHz ³¹P nmr spectroscopy of a cationic technetium complex for myocardial imaging

Dear Prof. Shapiro:

Technetium is the most desired isotope for nuclear medicine procedures. Until now there has been no Tc-99m agent available for determining myocardial perfusion function. Recently, we have prepared a cationic technetium complex that has been successfully used for myocardial imaging in a variety of animals including rabbits, dogs and baboons. The diamagnetic complex formulated $[Tc(dmpe)_3]^+$ has Tc in +1 oxidation state octahedrally co-ordinated by six phosphorus atoms; it is prepared by refluxing an aqueous ethanolic solution of NH_4TcO_4 with 5-molar excess of the ligand for 3 hours. The pale yellow complex is isolated as the chloride salt. Elemental analysis, FABMS and X-ray crystallographic studies confirm the structure(figure 1).

Following a report nearly two years ago on 99 Tc nmr of NH₄TcO₄(1), we were interested in looking at the 99 Tc and 31 P nmr of the above complex on our multinuclear Bruker WP-200. With a sensitivity of 0.376 relative to 1 H=1.0, we were able to establish quickly the observation parameters employing a D₂O solution of NH₄TcO₄. Because of the 30 KHz sweep-width limitation on our WP-200, we had to play with the synthesizer [setting=(200.133-X)/2] in order to observe the 99 Tc resonance of the Tc(dmpe)₃ complex(0.1 M in dmso-d₆); the septet (figure 2), resulting from coupling to six equivalent 31 P atoms, was observed in a short time. We back-calculated the chemical shift of the (center of the) 99 Tc resonance of the complex to be -1891 ppm w.r.t. NH₄ 99 TcO₄.

Observing the ^{31}P spectrum of the complex was straightforward; the ten-line pattern resulting from coupling to ^{99}Tc (I=9/2) and centered at 39.5 ppm w.r.t. external 85% H_3PO_4 was observed in a matter of minutes (Figure 3).

In all the above cases, proton noise-decoupling was employed to reduce linewidths and eliminate further splittings due to ¹H-X couplings which caused considerable broadening of the lines with fine structure at the top. With proton-decoupling, the line-widths at half-heights were about 25 Hz for the ⁹⁹Tc spectrum and 37.5 and 45 Hz, respectively, for the two outer lines and the eight inner lines of the ³¹P spectrum.

With best regards,

Butramanyay


V. Subramanyam

...

Karen Linder

Sincerely yours,

P. R. Srinivasan



Ref: 1. M. J. Buckingham, G. E. Hawkes and J. R. Thornback, <u>Inorganica</u> Chimica Acta, <u>56</u>, L41-L42(1981)



The Ohio State University

Department of Chemistry

140 West 18th Avenue Columbus, Ohio 43210 Phone 614 422-2251

March 31, 1983

Professor B. L. Shapiro
Department of Chemistry
Texas A and M University
College Station, Texas 77843

Alkynyllithium, rotational barriers in acyl dihydropyridines

Dear Barry:

In response to your dire notes I belatedly report the following news.

Alkynyllithium: We teach undergraduates that alkynylmetal compounds are salts containing alkynyl anions. The latter are stabilized thanks to the 50% plus "s" character associated with the sigma orbital on carbon. This idea, at least for lithium compounds is now thoroughly debunked by the results of an NMR experiment. Carbon-13 NMR of C_1 of 3,3-dimethyl-1-lithiobutyne in THF at -50° consists, for the ⁶Li enriched (96%) species, of a 7 line multiplet, closely, 1:3:6:7:6:3:1, and for the ⁷Li (normal) material 10 lines, 1:3:6:10:12:12:10:6:3:1 indicating three lithiums symmetrically coupled to each α -carbon. The coupling constants are $J(^{13}C, ^{7}Li)$ =16Hz and $J(^{13}C, ^{6}Li)$ =6Hz which is the expected ratio and confirms where the splittings are coming from. Altogether these results mean that this typical alkynyllithium compound exists as a cubic tetramer, where each lithium

is coordinated to one THF molecule. This is just what has been found for most simple organolithium etherates (methyl, vinyl, aryl). Further the ¹³C, lithium coupling constants are also very similar to the reported ones for other lithium compounds.

The alkynyllithium in THF is not any more ionic than the other RLi etherates and should not be described as a salt.

Dihydropyridines-barriers to rotation: Recently we became interested in using rotational barriers to assess the effect on bonding of inserting a heteroatom across the termini of a π -electron loop, viz 3-aza-1,4-dienyl, so we have barriers for the N-acetyl dihydropyridine, I, and its sulfone analog II. We think the

Professor B. L. Shapiro

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_4$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

$$CH_3$$

twisted state of II is stabilized (relative to twisted I) because the nitrogen lone pair and the double bonds we provide 6 π electrons <u>and</u> there is p,p overlap between carbon and sulfur, i.e. have an "aromatic" transition state for rotation.

I am interested in the comments of theoretically inclined readers.

Best wishes.

Yours sincerely,

Gideon Fraenkel

Professor of Chemistry

Post-Doctoral Position

There is now a Post-doctoral position open in my group for someone who will work on structure and dynamic behavior of organolithium compounds. This involves synthesis of isotopically enriched, ⁶Li and ¹³C, organolithium compounds and their study using ¹³C, ⁶Li and ¹H NMR line-shape analysis. A high field multinuclear instrument is used in this work as well as an electromagnet system assembled in this department. The person we appoint will also help in maintaining the electromagnet NMR.

Address inquiries to:

Professor Gideon Fraenkel Department of Chemistry Ohio State University 140 W. 18th Avenue Columbus, Ohio 43210 U.S.A.

GENERAL 🚳 ELECTRIC

CORPORATE RESEARCH AND DEVELOPMENT

GENERAL ELECTRIC COMPANY ● RESEARCH AND DEVELOPMENT CENTER ● P.O. BOX 8 ● SCHENECTADY, NEW YORK 12301 ● (518) 385-2211

April 4, 1983

Professor Barry L. Shapiro Department of Chemistry Texas A&M University College Station, Texas 77843

Carbon-13 NMR of PPO Copolymers

Dear Barry,

We have recently been examining some PPO copolymers in which the comonomer with 2,6-xylenol is 2-methyl-6-isopropylphenol. Variations in the polymerization conditions can produce polymers which range from fully random to highly blocked copolymers. Carbon-13 NMR is a very useful way of determining the structure of the final copolymer.

As shown in the expansion in the figure, the C-4 carbon of each ring resolves into two peaks as a result of sequencing. In the case of a random 50:50 copolymer (A=B) the spectrum shown is obtained in which each resonance is of the same intensity. In a more blocked copolymer, the intensities of BB and AA increase at the expense of the interblock units BA and AB.

We have examined polymers of different composition and with other alkyl groups in B and have found similar results. A full paper on these copolymers will be submitted soon.

Sincerely,

E. A. Williams

P. E. Donahue

Materials Characterization Operation R&D Applications Operation

/ldr




Figure. Carbon-13 nmr spectrum of PPO-Me₂/PPO-Me,i-Pr copolymer at 75.4MHz.

VYSOKÁ ŠKOLA CHEMICKO-TECHNOLOGICKÁ

Laborator synthetických paliv 166 28 PRAHA 6 DEJVICE, SUCHBATAROVA 5 Telefon 332/4168, 32 15 04 Dálnopis 122 744 - VSCH/C

Prof. Barry Shapiro
Department of Chemistry
Texas A&M University
College Station, Texas 77843


Praha, 6.4.1983 č.j. 1157/83

Title:

"Configuration and Conformation of 2-Alkyl-3, 3, 4-Trifluoroxetanes"

Dear Prof.Shapiro,

In 1976 we started together with Dr. Hrušovský who died accidently in summer of the 1980 year, a study of group of 2-alkyl-3,3,4-tri-fluoroxetames. This type of compounds is similar to the cyclobutane system previously studied by Ernst^{2,3} and Roberts^{4,5}.



$$R^1 = H, CH_3$$

 $R^2 = CH_3, C_2H_5, (CH_3)_3C, C_6H_5$

Now we have performed a detailed analysis of 2-alkyl-3,3,4-trifluor-oxetanes to determine the relative configuration of fluorine atom on carbon atoms 2 and 4 and to obtain the approximate conformation of the oxetane ring. The chemical shifts and the coupling constants were calculated from the experimental data by means of the iterative LAOCN3 procedure. Relative signs of the coupling constants were determined by INDOR and DQT experiments (the sign of $^2J_{CF}$ we supposed to be positive). From solvent effects, temperature dependences, LSR experiments (with Yb(FOD)₃) and a comparison of experimental data with the theoretical geometry of exetane ring calculated

by EHT method, we obtained following conclusions:

- 1) Conformation of oxetane ring is nonplanar.
- 2) Vicinal spin-spin coupling constant $^3J_{\scriptscriptstyle\mathrm{FF}}$ lies in the range from 8 to 12 Hz for trans configuration and from -4 to -11 Hz for cis configuration (of both fluorine atoms).
- 3) Fluorine atom in the CF2 group which is in pseudoequatorial position is deshielded in comparison with the pseudoaxial one.
- 4) If substituent R1=R2, the F atom in the dominant ring conformation is pseudoaxial. If $R^1 \neq R^2$, the conformation of the ring depends on the more bulky substituent, which occupies predominantly the pseudoequatorial position.
- 5) Fluorine atom F decreases the shielding of the gauche fluorine atom in the CF₂ group (probably due to an interaction with the electrical dipole of the CHF group).

We hope you will receive our contribution before writing us the second warning.

With best regards, yours sincerely

Petr Traka hil fleye (
Milan Hájok

Prague Institute of Chemical Technology

- 1. I.Hrušovský: Thesis, SVŠT Bratislava, 1977
- 2. R.R.Ernst: Mol. Phys. 16, 241(1969)
- 3. R.R.Ernst: J.Chem.Phys. 45,345(1966)
- 4. J.B.Lambert, J.D.Roberts: J. Amer. Chem. Soc. 87,3884(1965)
- 5. M. Takamashi, D.R. Davis, J.D. Roberts: J. Amer. Chem. Soc. 84,2935(1962)

UNIVERSITÄT GESAMTHOCHSCHULE SIEGEN

FACHBEREICH 8 Naturwissenschaften II (Chemie/Biologie) Organische Chemie II Prof. Dr. H. Günther

Universität-Gesamthochschule-Siegen · Postfach 101240 · 5900 Siegen

Prof. Bernard L. Shapiro Texas A&M University Department of Chemistry College Station Texas 77843 USA

5900 Siegen, April 6, 1983 Adolf-Reichwein-Straße Postfach 10 12 40 Telefon (02 71) 7401 Durchwahl 740-4390/4400

TANDEM-SEFT

Dear Barry,

having a lot of deuterated systems at hand — a consequence of our studies on $^{2}\text{H}/^{13}\text{C}$ isotope effects — we thought it would be fun to join the pulse game (1). Spin echo modulation by $^{13}\text{C},^{2}\text{H}$ coupling was a promising first try and the rules in one and two dimensions were quickly worked out (2). For practical applications, $^{2}\text{H-SEFT}$ is of interest for the analysis of deuterated compounds and we have developed, in combination with difference spectroscopy, a strategy for this purpose (3). The latest hit is, however, TANDEM-SEFT, where we pulse the ^{1}H and the ^{2}H decoupler independently:

13 C	π/ ₂ Γ τ	π	τ,	;	T2	TT I	T2	FID
¹H	7////			7			////	
² H	7////							1///

Minor instrumental modifications were necessary to perform this experiment on our WH 400 spectrometer. They are available on request. TANDEM-SEFT offers four unique possibilities for signal editing of partially deuterated systems not easily available otherwise:

- (1) Suppression of all proton- and/or deuterium-bearing carbons with sole detection of quarternary carbons [(τ_1 = 1/2J(C,H); τ_2 = 1/3J(C,D)];
- (2) Suppression of all deuterium-bound carbons [(τ_1 = 1/J(C,H); τ_2 = 1/3J(C,D)] and phase selection for CH_n groups;
- (3) Suppression of all proton-bound carbons [(τ_1 = 1/2J(C,H); τ_2 = 1/2J(C,D)] and phase selection for CD_n groups;
- (4) Combined selection of deuteron and proton multiplicity [(τ_1 = 1/J(C,H); τ_2 =1/2J(C,D)]. The details will be submitted to JACS.

Hope to see you in Chicago! Best regards,


D Cahmitt

T P Wesener

H. Günther

- (1) For a review see R. Benn and H. Günther, Angew. Chem., May or June '83.
- (2) P. Schmitt, J.R. Wesener, and H. Günther, J. Magn. Reson. 52, Nr. 3 (1983).
- (3) J.R. Wesener and H. Günther, Org. Magn. Reson., in press.

Introducing simplicity, versatility and productivity to NMR.



A new series of high resolution NMR spectrometers for even the most demanding experiments.

The unique features of the new AM Series provide truly new levels of simplicity and productivity along with the built-in Bruker versatility.

Check these against any other Supercon:

Advanced pulse program system, with software designed specifically for simplicity and clarity.

0 0

Complete microprocessor control of basic spectrometer functions, designed to simplify system set-up, lock operation and magnet shimming.

Fast 8-color display processor for display of complex 2-D

contour plots in seconds.

Fourier transform array processor, with associated

display and multi-pen plotter software for increased clarity and productivity.

- Large CPU memory capacity of up to 1 M word data, with 24 bit single precision arithmetic and a floating point instruction set.
- Extensive line of superconducting magnet systems from 200 to 500 MHz.
- Plus Bruker's flexible design concept, allowing for future upgrading.

So, if you want simplicity, versatility and productivity in your NMR work, write for complete information on the new AM Series or call a Bruker representative.

Bruker Instruments, Inc., Manning Park, Billerica, MA 01821. In Europe: Bruker GmbH, Silberstreifen,

D-7512 Rheinstetten 4, West Germany.


With the new Bruker AM Series.





JKER NMR Systems designed to solve problems.





In addition to the convenient weekly planner, the Bruker Almanac features over 30 tables of useful data, formulae and other information. It also contains a 50 page section on instrumentation: covering NMR, NMR Imaging, EPR, FT-IR, MS, LC, Magnets and Data Systems.

We have a copy reserved for you.

Almanac 1983

\square Please send me the Bruker Almanac 1983. FREI	E.
--	----

- ☐ I am interested in High Resolution NMR Systems (......MHz)
- ☐ Please send me information on the new AM Series described on the reverse.
- ☐ Better yet, please have your representative call me.

My Phone Number is: (_____) My Application is: ___

emo/Titlo:

Name/Title:

Institute/Company: ______

Department: ______

Address:_____

City/State/Zip:



Department of Chemistry

Colorado State University Fort Collins, Colorado 80523

April 11, 1983

Professor B.L. Shapiro Dept. of Chemistry Texas A&M University College Station, TX 77843

Note of Caution in Using Br-79 to Set the Magic Angle for C-13 MAS Experiments.

Dear Barry:

We have been using K 79 Br to set the magic angle for 13 C CP/MAS spectra for the past year (TAMU $\underline{279}$, 42, December 1981). The proximity of the resonance frequencies of 79 Br and 13 C is a great advantage in that the probe and power amplifier do not need to be retuned. In some 13 C gated-decoupled MAS spectra, however, we see a sharp signal that arises from 79 Br aliased into the 13 C region. At 3.5 T the resonances are only 140 kHz apart, and a 5 μs pulse excites both nuclides. The aliased 79 Br peak can be identified by opening the input filters or expanding the sweep width and noting the great increase in intensity.

The ¹H spin-temperature alternation feature of the cross polarization experiment effectively removes this FT signal, for we have not observed it even in long CP runs. Please credit this contribution to Gary Maciel's "account".

Sincerely,

James S. Frye

JSF:1b



LABORATORIUM VOOR FYSISCHE CHEMIE

Toernooiveld 6525 ED Nijmegen Telefoon (080) 55 88 33

FACULTEIT DER WISKUNDE EN NATUURWETENSCHAPPEN KATHOLIEKE UNIVERSITEIT NIJMEGEN, NEDERLAND

Prof. B.L. Shapiro Texas A & M University College of Science Department of Chemistry COLLEGE STATION, Texas 77843 U.S.A.

Uw kenmerk

Uw brief van

Ons kenmerk

U5259/AK,FS,WS/DvdW Datum

April 19, 1983

Onderwerp

Dear Professor Shapiro,

²⁷Al magic angle spinning NMR of zeolites

Crystalline alumino silicates like zeolites can absorb relatively small molecules like for instance water. Most of the water molecules are physisorbed in the pores of the material but some are chemisorbed at Si-O-Al and Si-OH sites. With H MAS NMR, water molecules at these two adsorption sites can be distinguished [1]. The Si-O-Al sites (Brønsted sites) play an important role in the catalytic activity of the zeolites and therefore the interaction between adsorbed molecules and these sites is of interest.

It is believed that water molecules adsorbed at Brønsted sites form clusters with the acidic proton: $\operatorname{H}^+(\operatorname{H}_2\operatorname{O})_n$ where n is of the order of two [2]. Therefore, hydration and dehydration of zeolites is expected to alter the position of the H^+ proton and consequently change the configuration of and charge distribution around the aluminum-oxygen tetrahedron. In a perfect tetrahedron the $^{27}\operatorname{Al}$ quadrupole interaction is zero and the $\frac{1}{2} \Leftrightarrow -\frac{1}{2}$ $^{27}\operatorname{Al}$ transition narrow. For a perturbed tetrahedron this NMR transition should broaden due to the quadrupole interaction, a broadening which is not eliminated by magic angle spinning [3]. A broadening by quadrupole interaction of the $\frac{1}{2} \Leftrightarrow -\frac{1}{2}$ transition is characterized by its field dependence, the larger the Zeeman interaction relative to the quadrupole interaction, the smaller the broadening.

Fig. 1 shows 27 Al MAS spectra of H-ZSM-5 at three different frequencies for the hydrated and dehydrated form. Clearly, on dehydration the 27 Al $_{2}^{1}$ \leftrightarrow $-_{2}^{1}$ NMR

LABORATORIUM VOOR FYSISCHE CHEMIE

Toernooiveld, Nijmegen Telefoon (080) 55 88 33

FACULTEIT DER WISKUNDE EN NATUURWETENSCHAPPEN KATHOLIEKE UNIVERSITEIT NIJMEGEN

Geadresseerde Prof. B.L. Shapiro Ons kenmerk U5259/AK,FS,WS/DvdW Datum 19-4-*83 Blad No

line broadens. The field dependence of the broadening establishes the quadrupole origin (at 15.6 MHz the line of the dehydrated sample is too broad to be detected). Clearly, removing the chemisorbed water changes the charge distribution around aluminum.

Fig. 2 shows ²⁷Al MAS spectra at 78.2 MHz for Na-ZSM-5 as a function of hydration. These spectra show not only a line from four-coordinated Al at 53 ppm, the genuine ZSM-5 aluminum, but in addition lines at 65 ppm and around 0 ppm. The lines around 0 ppm originate from six-coordinated Al (Al(H₂0)³⁺₆ in the pores and Al₂0₃ impurities). On dehydration the line from the four-coordinated Al of the catalytic active sites of Na-ZSM-5 broadens so much that it is practically unobservable. On rehydration it slowly reappears again. Apparently, the counter ion Na⁺ on dehydration perturbs the electric field gradient at Al more than H⁺ in H-ZSM-5. This makes it plausible that the increase of the electric field gradient at Al on dehydration is not only due to a change in the location of the counter ion H⁺ or Na⁺ but also due to a distortion of the aluminum-oxygen tetrahedra.


- [1] K.F.M.G.J. Scholle, W.S. Veeman, J.G. Post and J.H.C. van Hooff, Zeolites, in press.
- [2] H. Nakamoto and H. Takahashi, Zeolites, $\underline{2}$ (1982) 67.
- [3] E. Kundla, A. Samoson and E. Lippmaa, Chem. Phys. Letters, <u>83</u> (1981) 229.


Best regards,

A.P.M. Kentgens

K.F.M.G.J. Scholle

W.S. Veeman





770-10

Instruments Inc.

Orchard Park
PO Box 332 Danbury CT 06810

Telephone: 203 796 2500

April 7, 1983

Professor B. L. Shapiro Department of Chemistry Texas A & M University College Station, TX 77843

2D Magnetization Transfer in the Rotating Frame

Dear Barry:

We have for some time been interested in investigating chemical exchange in small peptides. The problem often encountered is that the resonances of interest are close together (<25 Hz) and the use of the more common techniques (transfer of saturation and selective inversion) becomes impractical. One instance of this is the cis-trans resonances of the C-4H histidine resonances in TRH (p-glu-his-pro-NH₂) where the separation for the two conformations is about 22 Hz (using an IBM Instruments WP270SY spectrometer).

Using the Hennig and Limbach technique of magnetization transfer in the rotating frame (J. Magn. Resn. 49, 322-28 (1982)) we measured the $T_{1\rho}$ of the two resonances using the fixed delays $\tau = \frac{1}{\Delta f}$ and $\tau = \frac{1}{2\Delta f}$ with a frequency offset of Δf (22.06 Hz) downfield from the trans resonance. Since spin-locking effectively narrows the lines in the rotating frame the technique is potentially useful in cases where the resonances under investigation have only a small chemical shift difference.

After completing the one dimensional experiments we decided to extend the experiment to two dimensions (J. Jeener et. al., J. Chem. Phys. 71(1), 4546 (1979)). The fixed delay τ in the one dimensional experiment becomes the incremented delay DO (and the F1 frequency domain after Fourier transformation). The initial value of DO was 1 msec and was incremented by 2.5 msec 63 times. The 90° (P1) pulse was 9 μ sec and the spin locking (P2) pulse was 200 msec. This pulse (P2) corresponds to the mixing time in 2DNOE experiments. The two pulses had a relative phase shift of 90°. Sixteen transients were acquired in each of the 64 files.

The contour plot shown opposite reveals the results of the experiment. The 'wings' at the bottom of the plot are present because constant phase was used in the F1 dimension and only the region from 0 to -90 Hz is plotted. We performed the 2DNOE experiment on this molecule at 60°C and saw cross peaks. However, for the two dimensional Hennig and Limbach experiment it was necessary to raise the temperature to 75°C to see the cross peaks. This is probably due to the fact that the T_{10} is shortened by coupling to N-14 which is modulated by quadrupolar relaxation (Bleich and Glasel, Biopolymers, 17, 2445 (1978)).


Sincerely,

Joyce A. Wilde

Hermann E. Bleich

Richard E. Schwemm

/lw Enc.



University of East Anglia

Professor R.K. Harris

School of Chemical Sciences University of East Anglia Norwich NR4 7TJ, England

Telephone Norwich (0603) 56161 Telegrams UEANOR NORWICH

April 13 1983

Professor B.L. Shapiro Dept. Chemistry Texas A&M University College Station,TX 77843 U.S.A.

Dear Barry,

2-D Silicon-29 Spectra

Your green letter arrived on my desk at the same time as some splendid 2-D spectra obtained using the Bruker 400 MHz spectrometer (ca.80 MHz for ²⁹Si) at the University of Warwick, so we are able to respond to your prodding quickly. The chemical problem, which has been fascinating us for some time, is that of establishing the structures of the complex mixture of anionic species which are in equilibrium in alkaline silicate solutions. made extensive use 1,2 of enrichment with 29Si, in combination with tedious point-by-point homonuclear decoupling experiments. The full spectra are complicated, so 2-D work offers elegant alternatives. The attached partial J-resolved spectrum shows the power of these techniques. We also have useful COSY spectra of these systems. We believe these are the first 2-D homonuclear ²⁹Si spectra in existence. The notation for the structures given on the figures is such that Q^1 , Q^2 and Q^3 represent silicate units with one, two and three siloxane bridges respectively. So far, unfortunately, although we've nicely confirmed some previous deductions, no really new information has emerged - but we're still hoping.

Best wishes,

Yours sincerely,

Kovin


ROBIN HARRIS

MAX O'CONNOR*

- 1. R.K. Harris, C.T.G. Knight & W. Hull, J. Amer. Chem. Soc. 103, 1577 (1981)
- 2. R.K. Harris, C.T.G. Knight & W.Hull, ACS Symp. 194, 79 (1982)

^{*} On leave from La Trobe University, Melbourne.

PARTIAL 2-D J-RESOLVED ²⁹Si NMR SPECTRUM of an alkaline potassium silicate solution ~1.41 M in SiO₂ K:Si = 1:1 278 K



LABORATOIRES DE CHIMIE

Département de Recherche Fondamentale de Grenoble

N/référence DRF/CH/83-134/AP

Grenoble, le April 5, 1983

Y/réf.: BLS/lmk March 9th, 1983

Professor B.L. Shapiro, Department of Chemistry, Texas A M University, College Station, TEXAS. 77843. U.S.A.

NMR Investigation of the cellulose conformation in an aprotic solvent solution

Dear Professor Shapiro,

The recent discovery of new organic solvents (1) for cellulose, allows to investigate the conformation of this natural macromolecule in solution, by using ¹H NMR. Among these solvents NN-dimethylacetamide containing 5 to 10 % lithium chloride, is of particular interest, due to its aprotic character. As shown by Figure 1, it is possible to observe the hydroxyl proton signals of cellulose in this solvent. The low signal resolution of this spectrum is due to the high viscosity of the solution.

Figure 2 shows the temperature effect on the hydroxyl proton signals. The chemical shift increment with temperature are:

$\frac{\Delta\delta}{\Delta T}$ 10 ⁻³ ppm/°C	1.9	2.5	1.	6.75	
Proton	H-1	OH-2	OH-3	OH-6	

It is known that the effect of temperature on proton MR signals is mainly governed by a change in intermolecular and intramolecular association. The lower value of $\Delta\delta/\Delta T$ for OH-3 observed here has also been found for cellobiose, the cellulose dimer, in a serie of aprotic solvents (2). Generally, hydroxyl group involved in intramolecular hydrogen found are less affected by temperature variation. To explain the particular behaviour of OH-3 we propose an intramolecular hydrogen bond between OH-3 and the ring oxygen atom 0-5' of the next D glucopyranosyle unit of the chain, as it has been found in the solid state.

Sincerely yours,

Dr Marc VINCENDON

- into




Figure 1 : ¹H NMR spectrum (250 MHz) of cellulose in N,N-dimethylacetamide-LiCl solution at 60°C.




Figure 2: Temperature variation of cellulose hydroxyl protons with temperature.

References:

- 1) A.F. TURBAK, A. EL KAFRAWY, F.W. SNYDER and A.B. AUERBACH French Patent 2, 466, 473 (1981)
- 2) D. GAGNAIRE, J. SAINT GERMAIN and M. VINCENDON J. Appl. Polymer Sci in Press
- P.S.: Please credit this contribution to Dr Roselyne BRIERE account.



Sundance Supply Company

P.O. Box 137 • Sunol, CA 94586 • (415) 862-2021

EQUIPMENT FOR SALE

Spectrometer NMR XL-100-15 system with: pulse preamp, controller, radio amplifier, data machine console, TTY, wide band amplifier, storage scope, and RF power amp.

12" and 15" Varian NMR Electromagnets, Maximum rating 14 and 23.5 Kilogauss.

For more information, call Clete Hillseth at (415) 862-2021.



University of Illinois at Chicago

HEALTH SCIENCES CENTER

Graduate College
RESEARCH RESOURCES CENTER
Bioinstrumentation, Biostatistics, Electron Microscope,
Environmental Stress, Instrument Shop, Scientific
Computer, and Spectroscopy Facilities

1940 West Taylor Street Chicago, Illinois 60612 (312) 996–7600 Mailing Address: P.O. Box 6998 Chicago, Illinois 60680

April 6, 1983

Dr. Bernard L. Shapiro Department of Chemistry Texas A & M University College Station, TX 77843

Dear Dr. Shapiro:

Recent months have been exciting. The University of Illinois at the Medical Center (UIMC) and the University of Illinois Circle Campus (UICC) have been consolidated since September 1, 1982 into one entity called the University of Illinois at Chicago. This new institution has two campuses referred to as the UIC Health Sciences Center (HSC) and UIC University Center (UC), respectively.

The Research Resources Center of the HSC has recently purchased two (2) Nicolet broadbanded NMR spectrometers—a 360 NB and a 200 WB. These instruments are positioned in the basement of our building to avoid the vibration problems experienced by some instruments on upper floors of a building. Of course, during the installation it was discovered that we had a floor vibration problem. Vibrations above 20 Hz were not too troublesome but vibrations below 10 Hz made shimming for line shape very difficult. The troublesome side bands looked like spinning side bands except that they would vary in frequency and were not related to spinning rate.

To prove that these artifacts were caused by just floor vibrations, we measured just floor vibrations using a blood pressure transducer with its diaphragm resting on a metal washer set on the base plate of the magnet. Only floor vibrations caused the diaphragm to oscillate which was observed using a carrier preamplifier with an oscilloscope or strip chart recorder. The frequency offset of the troublesome side bands correlated very well with the frequency of the floor vibrations.

To reduce the problem, the energized magnet was lifted off its base plate, four (4) rubber air cushions were installed under the base plate, and the magnet was lowered and reconnected to the base plate. The total cost in dollars was \$350. The total cost in aggrevation is not printable. The problem was reduced by about a factor of 100.

Michael Mutaw

Biomedical Engineer

Sincerely

Robert A. Kleps, Ph.D. Senior Spectroscopist

RAK:ac

PHYSICAL CHEMISTRY LABORATORY OXFORD UNIVERSITY

Telephone OXFORD (0865–) 53322



SOUTH PARKS ROAD OXFORD OX1 3OZ

31st March, 1983

RF/ml

Professor Barry Shapiro, Department of Chemistry, Texas A and M University, College Station, Texas 77843, U.S.A.

Dear Barry,

Practical tips for Broadband Decoupling

Most readers will be aware of the recent new methods of heteronuclear broadband decoupling using "supercycles" (1) built up from magic cycles (2) of composite spin inversion pulses (3). Some results were reported in this Newsletter about a year ago (TAMUNMR 284 May, 1982). The purpose of this note is to point out some of the experimental precautions needed to obtain the very best decoupling performance and to advocate a more recent scheme called WALTZ-16 (4,5). We emphasize that for many routine carbon-13 applications (e.g. where very high resolution is not important) these adjustments may well be superfluous.

1. Inaccurate radiofrequency phase shifts.

When observing carbon-13 line widths of the order of 0.1 - 0.2 Hz, sequences which use 90° radiofrequency phase shifts are quite sensitive to the accuracy of the three phase shifts involved. A simple solution to this problem is to use the spin inversion element (4)

$$R = 90^{\circ}(x) 180^{\circ}(-x) 270^{\circ}(+x) \equiv 1 \overline{2} 3$$

All sequences which employ only 180° radiofrequency phase shifts can be shown to be insensitive to errors in the phase shift (5). This element is therefore superior to earlier composite inversion pulses, simpler to implement, and covers a rather wider proton bandwidth. We have gone over completely to the WALTZ-16 sequence. With $\gamma B_2/2\pi=2 \text{kHz}$ (about 0.5 watts of power supplied to the Varian XL-200 probe) it covers a proton range of approximately 4 kHz. Incidentally, at these low power levels, we have noticed that shimming the B field is a much easier operation.

2. Amplitude Imbalance between Channels

With the WALTZ-16 scheme there are only two decoupler channels (0° and 180°). Computer simulations and experiments suggest that the B₂ level in each channel should be equal within better than 0.5% for optimum decoupling. If the phase shift is generated at an early stage of the amplifier chain (as in Varian XL-200) good amplitude balance is achieved automatically. However if the decoupler phase shift is accomplished by way of a double-balanced mixer,

there could be difficulties, either because the balance is insufficiently precise or because the rise and fall times of the driving pulse amplifiers are too slow. The cure is to use a current amplifier (for example National Semiconductors LH0002C) to drive the double-balanced mixer. No decoupling capacitor is used on the output since this would slow down the switching. A suitable circuit is shown in Figure 1. Setting up the amplitude balance is first done roughly by observation on an oscilloscope. Then the decoupled peak height is monitored for a very narrow line (few tenths of 1 Hz) for an offset $\Delta B = 0.25$ B₂ (chosen because it is relatively sensitive), carefully adjusting the 10-turn balance potentiometer for optimum peak height. This is a once-and-for all adjustment.

3. Cycling Sidebands

The cycling rate of the decoupler is much slower than the carbon-13 sampling rate, and this introduces very weak "cycling sidebands" into the decoupled spectrum. Many of these have random phases if the proton pulsing and carbon sampling processes are desynchronized. The cycling sidebands can be considerably attenuated by using a "blackbox" hardware decoupler instead of computer-controlled decoupler pulsing. A suitable circuit (5) is shown in Figure 2.

Kindest regards,

James Keeler A. J. Shaka Ray Freeman

Ray

^{1.} Levitt, Freeman and Frenkiel, J. Mag. Reson., 50, 157 (1982).

^{2.} Levitt and Freeman, J. Magn. Reson. 43, 502 (1981).

^{3.} Levitt and Freeman, J. Magn. Reson. $\overline{33}$, 473 (1979).

Shaka, Keeler, Frenkiel and Freeman, J. Magn. Reson. 52, 335 (1983).

Shaka, Keeler and Freeman, J. Magn. Reson., to be published.

Figure 1






Figure 2



Only one NMR instrument surpasses the sensitivity and performance of Varian's XL-300...

The new XL-400

featuring
Varian's new
high-speed,
color-graphic
ADVANCE™
NMR Data
System



Varian's new ADVANCE XL-400 Superconducting FT NMR Spectrometer gives you:

Guaranteed 400:1 S/N. That's the highest specification in the industry, on ¹³C ASTM 10mm, regardless of magnetic field strength.

8K floating point transforms in less than 500 milliseconds. When equipped with its optional array processor, our new ADVANCE Data System performs Fourier transforms on the XL-400 faster than any competitive instrument. Results are far more accurate, too, with Varian's 32-bit acquisition processor and floating point mathematics.

See more information—faster—on the 16-color CRT. The ADVANCE XL-400's operator interface, with its large display screen and wide spectrum of colors, allows you to make quicker and more accurate decisions based on CRT visuals.

More memory than any other NMR system.


The new ADVANCE Data System features a multicomputer design with 464K of memory and dual processors, each further expandable to 16 megabytes: big enough to handle the largest 2D data tables.

Flexible and easily expandable. This system lets you add commercially available peripherals to meet specific application requirements. This capability and our expandable Pascal-based software keeps your instrument "state-of-the-art" in the rapidly evolving field of NMR research.

Send for literature now. For details concerning Varian's new ADVANCE XL-400 Spectrometer, call the Varian sales office nearest you. Or write: ADVANCE XL-400, Varian Associates, D-070, 611 Hansen Way, Palo Alto, CA 94303.

A new ADVANCE in NMR ...from VARIAN

The ADVANCE XL-400 delivers these exciting results...




MM X20

Write or call today
to receive new
Illerature on the
ADVANCE NMR
Data System and
the XL-400
Spectrometer.

whymph X10

with NOE 300:1







Papanicolaou Cancer Research Institute at Miami, Inc.

1155 Northwest 14th Street • Miami, Florida • Telephone (305) 324-5572

Mailing Address: Box 016188 • Miami, Florida 33101

Prof. Bernard L. Shapiro Department of Chemistry Texas A & M University College Station, Texas 77843 April 14, 1983

Dear Prof. Shapiro:

In-Vivo ³¹P NMR Spectrum of a Rat Tumor


I wish to report that we have entered a phase of our work which involves obtaining ^{31}P NMR spectra from tumors in live animals. The spectrum shown here was obtained at 36.2 MHz from a rat mammary adenocarcinoma using a 10 mm coil on the tumor and 1600 54° pulses. Convolution difference was carried out with 15 Hz and 400 Hz line broadening factors. Preliminary spectral assignments are shown where SP = sugar phosphates, P_i = inorganic phosphate, PCr = phosphocreatine, NAD = reduced nicotinamide adenine dinucleotide, ATP and ADP = adenosine tri- and di-phosphate respectively.

I would appreciate it if you would reinstate our subscription to the Texas A & M NMR Newsletter. Thank you very much.

Sincerely yours,

Ronald E. Block, Ph.D. Associate Scientist

Ros Block



StonyBrook

Professor Bernard L. Shapiro Department of Chemistry Texas A & M University College Station, Texas 77843 Department of Chemistry
State University of New York at Stony Brook
Stony Brook, New York 11794
telephone: (516) 246-5050

April 14, 1983

FACTS ABOUT AQUEOUS SHIFT REAGENTS FOR HIGH RESOLUTION CATIONIC NMR

Dear Barry:

We have continued our work (J. Mag. Res., 46, 348 (1982)) on the above named shift reagents (SR). The accompanying Figure shows the isotropic hyperfine shifts, Δ , induced in several nuclei by $\text{In}(\text{CA})_3^{6-}$ SR (CA³⁻ is chelidamate, 4-hydroxypyridine-2, 6-dicarboxylate). The Δ values were calculated as the differences in observed reasonant frequencies in the presence and absence of SR. The spectrometers were locked on the $^{2}\text{H}_2\text{O}$ resonance of the solvent.

The ratios of the shifts induced in the 23 Na⁺ resonance by Dy(CA)₃⁶⁻ over those induced by Tm(CA)₃⁶⁻ (bottom two curves) are quite close to the value of -1.9 expected for a pure dipolar (pseudocontact) mechanism. This is strong evidence that this mechanism dominates, as expected, and that field-frequency locking eliminates most of the bulk-susceptibility contributions. We have observed similar results with the Dy(TTHA)³⁻/Tm(TTHA)³⁻ (TTHA⁶⁻ is triethylenetetraminehexaacetate) and Dy(PPP)₂⁷⁻/Tm(PPP)₂⁷⁻ (PPP⁵⁻ is tripolyphosphate) pairs.

The facts that the $^{25}\text{Mg}^{2+}$ and $^{39}\text{K}^+$ shifts induced by Dy (CA) $^{63-}$ are significantly greater than that of $^{23}\text{Na}^+$ mean that the binding of the former ions is greater and/or that the spatial coordinates of the binding sites are different. The ^{23}Na curve can be fitted quite well with the 1:1 thermodynamic binding constants $K_{\text{Na}} = 1579~\text{M}^{-1}$ and $K_{\text{Li}} = 61.5~\text{M}^{-1}$ and a limiting shift $^{60}\text{Na} = +11.9~\text{ppm}$. The fact that observed shifts of $^{25}\text{Mg}^{2+}$ and $^{39}\text{K}^+$ rise significantly above +12 ppm is strong evidence that the binding sites are different. (The field-independence of some ^{23}Na shifts has been established.) Limiting dipolar shifts (in ppm) for identical binding sites would be independent of the nucleus observed. We observe similar behavior for the $^{39}\text{K}^+/^{23}\text{Na}^+$ pair with Dy (TTHA) and Tm (TTHA) $^{3-}$. We find an extreme example of binding site difference in that Dy (NTA) $^{2^{3-}}$ (NTA) is nitrilotriacetate) shifts the $^{43}\text{Ca}^{2+}$ resonance upfield but the $^{113}\text{Cd}^{2+}$ resonance downfield. There are, of course, also differences in binding. The binding constants for Na⁺ and Li⁺ to Dy (CA) $^{3^{-}}$ reported above show this. We have also found that K⁺ out-competes Na⁺ for binding to Dy (PPP) $^{2^{-}}$ (the Guptas' SR).

In these kinds of experiments, there are three essential components: a) the In ion, b) the metal cation observed, and c) the chelate ligand. Examples of individual variation of a) and b) are given above. As an example of the variation of c), we report that while Dy(CA)₃6- shifts the ²³Na⁺ resonance upfield (Figure), Dy(TTHA)³⁻ shifts it downfield. More details of much of this work will be appearing soon (July) in <u>Inorganic Chemistry</u>.

Please credit this contribution to the Stony Brook subscription (c/o Paul Lauterbur).

Best regards,


Charles S. Springer, Jr. Assoc. Prof. of Chem. James a Balschi Research Assist. David M. Yarmush Postdoctoral Assoc.

Dave Garmush

Simon C. Chu Research Assist.

Mass. Inst. of Tech.

(Also contributing to this work were Marty Pike, Dave Hillman, Bob Lenkinski, Jean DeLayre, Eric Fossel, and Tom Smith).



The isotropic hyperfine shift, Δ , of various metal cationic resonances versus the stoichiometric mol ratio of shift reagent to cation, ρ. all cases, the metal chloride (MgCl2, KCl, NaCl) concentration is held constant at 50 mM. The magnetic isotope observed is present at natural abundance. In all cases, tris (tris(hydroxymethyl)aminomethane) buffer was present at a total tris concentration of 10 mM. The dashed curves are intended merely to guide the eye. For the upper three curves, the shift reagent was Li₆Dy(CA)₃·3LiCl. In the MgCl₂ case (²⁵Mg NMR, 24.5 MHz, 9.40T), the pH varied smoothly from 8.4 in the solution with no shift reagent to 7.8 in the solution with 150 mM SR. The temperature was ca.297°K. In the KCl case (39 K NMR, 14.0 MHz, 7.05T), the pH varied from 8.7 to 7.7. The temperature was ca. 297°K. In the NaCl case (23 Na NMR, 26.5 MHz, 2.35T), the pH varied from 8.6 to 7.9. The temperature was ca. 301°K. For the lower curve, (23Na NMR, 26.5 MHz, 2.35T), the shift reagent was Li6Tm(CA)3.3LiCl. The pH varied from 8.9 to 7.9. The temperature was ca. 301°K.



National Institutes of Health - Bethesda, Maryland 20205 \

April 15, 1983

Dr. B. Shapiro
NMR Newsletter
Department of Chemistry
Texas A&M University
College Station, Texas 77843

Title: NMR studies of perfused cells; Mag. Res. Biol. Vol. 2.

Dear Barry:

I would like to briefly describe a novel technique which we have developed for the perfusion of cells. The advantages of our technique are that it is applicable to any cells grown in culture and can be used with any standard NMR spectrometer. The method involves embedding the cells in an agarose gel matrix 1 extruded in the form of a fine thread. This allows rapid diffusion (ca. 1 min) of metabolites, or drugs, into the cells. We are not aware of any other method which has these advantages. We have published a preliminary description of the method 2 and applications to mammalian cells using ^{31}p and ^{1}H NMR are in progress (Figure). 3

I also take this opportunity to inform you about the publication of Vol. 2 of Magnetic Resonance in Biology. This contains excellent chapters by: Debrosse & Villafranca, on "Isotope effects on phosphorus chemical shifts: applications to enzyme mechanisms;" Makinen & Kuo, on "Spin label probes of enzyme action" and Shindo, on "DNA backbone conformation and dynamics." We have also written a chapter (Lou Hughes, Jan Wooten and myself), entitled "Observations of amino acid side chains in proteins by NMR methods," which focusses on individual side chain types (His, Met, Tyr, etc.) and contains ca. 350 references. For those who are interested in water in tissue, and its potential relevance to NMR imaging there is a chapter by Lynch, on "Water relaxation in heterogeneous and biological systems." Altogether a varied yet in-depth collection. Additionally I draw your attention to two other reviews on different topics that we have written this year (a busy one), namely on "31p NMR studies of DNA and RNA", and a general review entitled "NMR in Biology and Medicine" being published as a single tract.

In further news from our laboratory I should emphasize we are still very active in ³¹P NMR studies of DNA conformation, particularly in relation to the B to Z transition. And lastly I confirm that I will be transferring to the National Cancer Institute and setting up a new laboratory there devoted to the multivariate uses of our favorite technique.

Yours sincerely,

Jack S. Cohen

Laboratory of Theoretical and


Physical Biology

National Institute of Child Health and Human Development

References

- 1. L. Jacobson and J. S. Cohen Biosci. Reports 1 141 (1981).
- 2. D. Foxall and J. S. Cohen J. Mag. Res. 51 346 (1983).
- 3. R. Knop, C. Chen, D. Foxall, J. S. Cohen, J. Mitchell and A. Russo, submitted.
- 4. J. S. Cohen (ed) Magnetic Resonance in Biology, Vol. 2, Wiley-Interscience, 1983.
- 5. C. Chen and J. S. Cohen, to appear in "Phosphorus-31 NMR: Principles and Applications," D. Gorenstein (ed), Academic Press, in press.
- 6. J. S. Cohen, R. Knop, G. Navon and D. Foxall, <u>Life Chemistry Reports</u>, Harwood Academic Press, in press.
- 7. C. Chen, J. S. Cohen and M. Behe, Biochemistry, in press.

Figure Legend: A stacked plot of 31 P NMR spectra at 109 MHz of Chinese hamster lung fibroblasts (ca. 108 cells) perfused in agarose gel threads at a flow rate of 2.7 ml/min with F12 medium (plus 10% FCS, 20mM Hepes, pH 7.35) in a Wilmad 10 mm screw cap tube. Each accumulation took 22 min (77° pulse, 220 msec delay) at 25°C. On stopping perfusion, the ATP signals decrease while the P_i signal concomitantly increases. On reperfusing the reverse occurs and the original levels are attained. This indicates that the cells can be kept metabolically active under the conditions of the experiment over a period of many hours. Microscopic examination, plating and other criteria indicate the cells remain morphologically normal.



296-35

DIVISION OF PHYSICAL CHEMISTRY THE ROYAL INSTITUTE OF TECHNOLOGY

S-100 44 STOCKHOLM 70

SWEDEN

Cable address: Technology

Stockholm, April 19, 1983

Professor B L Shapiro Texas A&M University College Station TX 77843 USA

Dear Professor Shapiro,

SIMULTANEOUS MEASUREMENTS OF SELF-DIFFUSION COEFFICIENTS AND SPIN-LATTICE RELAXATION TIMES

Measurements of self diffusion coefficients with the FT-NMR pulsed-gradient spin-echo method (1) have proved to be a valuable tool for elucidating structures of complex liquid systems, e.g. microemulsions and micellar solutions (2). The pulse sequence used is shown in Fig. 1a.

By combining the inversion recovery-method for measuring T₁ with the pulsed gradient spin-echo method, as shown in Fig. 1b, we have found it possible to simultaneously measure T_1 as well as the self-diffusion coefficients. The attenuation is in this case given by

$$\frac{A}{A_0} = (1 - be^{-\tau/T} 1) e^{-\gamma^2 g^2 D \delta^2 (\Delta - 1/3\delta)}$$
 [1]

where A_0 is the amplitude of a line in the Fourier transformed spin-echo without field gradient and A the corresponding amplitude with gradient pulses present, g is the magnitude of the field-gradient, D is the self-diffusion coefficient and τ , Δ , δ are time-parameters defined in Fig. 1b.

Fig. 2 shows the results from a typical ¹H measurement. The sample used was an un-degasser 1:1 by weight mixture of C6H6 (Merck, p.a.) and (CH3)2CO (Merck, p.a.).

The measurements were performed at 27°C with a Bruker CXP-100 spectrometer operating at 90 MHz, using quadrupolar gradient coils (3). The field gradient was determined by calibration with a H2O-sample; the calibration yielded a value of g = 1.30 ± 0.01 G/cm. T₁ and D were evaluated by a nonlinear curve-fitting of the parameters in Eq. [1] to the data. The values are consistent with those obtained in separate measurements, with comparable accuracy for the self-diffusion coefficients and somewhat lower for T_1 . The idea behind this method can, of course, also be applied to simultaneous measurements of T_1 and T_2 by combining the "inversion recovery" experiments with a CPMG-sequence.

Sincerely yours, Torbjörn Wärnheim

References:

P Stilbs and M E Moseley, Chem. Scr. 15, 176 (1980). B Lindman, P Stilbs and M E Moseley, J. Colloid Interface Sci. 83, 569 (1981).

3.a H Breuer, Rev. Sci. Instr. 36, 1666 (1965);

b G ödberg and L ödberg, J. Magn. Reson. 16, 342 (1974).




Fig. 1 a) Pulse sequence used in the pulsed-gradient spin-echo experiment.

b) Pulse sequence for simultaneous measurements of self-diffusion coefficients and spin-lattice relaxation times.




Fig.2 Spectra for 1:1 benzene-acetone mixture at different values of τ and δ . The time interval Δ =0.140 s in all experiments.

296-37STATE UNIVERSITY LEIDEN - THE NETHERLANDS GORLAEUS LABORATORIES - DEPARTMENT OF CHEMISTRY

P.O. Box 9502, 2300 RA Leiden

Phone: 71-148333, extension:3921

Your letter:

Your ref.:

Our reference:

Subject:

Dr. J. Lugtenburg

LEIDEN. April 20, 1983.

Prof.dr. Bernard L. Shapiro
Texas A&M University
College of Science
Department of Chemistry
College Station, Texas 77843
U.S.A.

Dear Professor Shapiro,

 1 H- and 13 C-NMR spectra of 11 13 C all-trans-retinal (1). For our studies in the rhodopsin and bacteriorhodopsin field we prepared all-trans-retinal with 92% 13 C incorporation on carbon position 11. The source of 13 C label was 1 13 C acetonitrile (M.S.D. 92%).




fig. 1 $11-\frac{13}{11}$ C-all-trans-retinal (1)

In figure 2(A) the 300 MHz ¹H-NMR spectrum of 1 is depicted.




fig. 2 300 MHz $^{1}\mathrm{H-NMR}$ spectrum of 1

In 2(B) the vinylic region is expanded and in 2(C) the same region of all-trans--retinal is recorded. From the spectrum the value for $J_{C_{11}-H_{11}}$ =150.4 Hz is obtained and the signals at 7.14 ppm show the presence of 8% unlabeled and 92% 13 C labeling at C $_{11}$. All other signals are essentially unchanged. 1 In figure 3A the 1 H noise decoupled 13 C-NMR(75.5 MHz) spectrum is shown. 2 In 3B the strong singlet of 1 at 132.5 ppm corresponding to the 92% incorporation at C_{11} is clearly visable.




fig. 3

From the natural abundance $^{13}\mathrm{C}$ NMR spectrum of 1 the following $^{13}\mathrm{C}$ - $^{13}\mathrm{C}$ coupling constants could be obtained.

 $\begin{array}{ccc} c_9 - c_{11} = 0 & c_{11} - c_{14} = 8.2 \text{Hz} \\ c_{11} - c_{13} = 0 & c_{11} - c_{19} = 4.3 \text{Hz} \end{array}$ $C_{10} - C_{11} = 58.7 \text{ Hz}$

C₁₁-C₂₀=3.3Hz

 $C_{11}-C_{12}=69.8 \text{ Hz}$

Up to now in the hydrocarbon field only the $^{1}\mathrm{J}$ $^{13}\mathrm{C}$ - $^{13}\mathrm{C}$ values of a few simple molecules are known. This is the first case where a larger polyene chain has been studied. The values for the various types of bonds agree very well with those known for small systems. We hope to prepare via 1 rhodopsin and bacteriorhodopsin with 92%. ¹³C on position 11 in their chromophoric part. Magic Angle Spinning ¹³C NMR spectroscopy of these molecules is expected to give novel essential information about the interaction of chromophere and apoprotein in these systems.

Sincerely yours,

Patrick

(J.A. Pardoen)

(C. Erkelens)

1) For the 1 H-NMR values of all-trans-retinal, D.J. Patel, Nature $\underline{221}$, 825-828 (1969). 2) For the 13 C-NMR values of all-trans-retinal, G. Engbert, Helv. Chim. Acta, $\underline{58}$,

2367-2390 (1975).

University of Bristol

School of Chemistry

19th April1983

Cantock's Close Bristol England BS8 ITS

Telephone: Bristol 24161 Ext.

Dear Professor Shapiro,

We do a lot of work with diphosphorus compounds, so come across AA'X spectra (A = 31 P, X = 13 C) quite often. They are also common in coupled 13 C spectra where A = 1 H.

The typical symmetrical 5-line spectrum (see figure) is often simplified to 4 lines (if $_{AA'}$ << L, $_{Imid} \approx 0$) or 3 lines (if $_{JAA'} >>$ L, $_{Iout} \approx 0$). The first case occurs frequently when A = $_{I}^{I}$ H and the carbon is bonded to one of the protons; the second is common when A = $_{I}^{31}$ P.

When the 5-line pattern is observed it is obviously impossible to extract the three coupling constants from the two line separations involved, although in many P-P cases one of the coupling constants can be assumed to be zero. We have therefore investigated the use of line intensities to provide further information. The following simple relationships ar used:

$$J_{AA'} = 2D \sqrt{I_{mid}/2} = 2D \sqrt{1 - I_{out}}$$

This technique has been tested for a number of cases where $A = {}^{1}H$, and compared with the results obtained by the use of ${}^{13}C$ satellites (see Table). The carbons that give five-line spectra are not those bonded to the two protons, but those for which $L \approx J_{AA}$.

For diethylmethylenemalonate the satellite method is unusable, but the 5-line spectrum of the carbonyl carbons (see Figure) gives a value of 0.88 ± 0.05 Hz for the geminal H-H coupling constant. The alternative method of determining this coupling constant, deuterium exchange, would probably fall with such small couplings, quite apart from the synthetic difficulties of preparing the labelled compound.

A problem that can arise when J_{AA} , < L is that the N doublet and the combination lines are of such similar intensity that they cannot be assigned. In this case the off-resonance technique of Fritz and Sauter (J. Mag. Res. 1974, 15, 177) can be used.

Please credit this contribution to Robin Goodfellow's subscription.

Yours sincerely.

Martin Murray

Kevin Higgins
Kevin Higgins

Figure: ¹³C Spectrum of Diethylmethylenemalonate




Table: Values of J_{H-H} for Symmetrical Molecules

from 13 C of \underline{C} from 13 C satellites of \underline{H}

11.8 + 0.2

$$\frac{H}{4}C - C = C + \frac{H}{3}$$
 3.39 ± 0.07
 3.5 ± 0.2
 $\frac{H}{0 - C}C = 0$
 5.74 ± 0.07
 5.7 ± 0.1
 $\frac{H}{000}C = C = 0$
 $\frac{H}{000}C = 0$
 $\frac{H}{00$

11.74 <u>+</u> 0.18



Department of Chemistry

4/8/83

Fort Worth, Texas 76129 817-921-7195

Dr. B.L. Shapiro Department of Chemistry Texas A&M University College Station, TX 77843

Dear Barry:

It is now clear that twenty years of contributions to the Newsletter does not absolve one from your dunning letters. I am thus reassured that nine months after I ultimately depart this earth someone will remember me.

We have been studying the electrophilic additions of several types of reagents to the benzobicyclol2.2.2 octadienes below for the past several months. Depending on the nature of the reagent and of the R groups one can get either syn or

anti addition of the reagent in either a syn or anti sense with reference to the aromatic ring side of the double bond. Often one gets rearrangements. Given the complexities of the situation, the use of NMR for structure assignment would be helpful. Even though carbon assignments can often be made by conventional means, it is not always evident whether rearrangement has occurred or not.

Aromatic methoxy or acetoxy groups bind Yb(FOD)₃ very weakly when compared with hydroxy, epoxy of even keto groups on the bicyclic ring. An exception is the rearranged iodohydrin below where the geometry allows the molecule to serve as a strong bidentate ligand. The values in () are the relative It would be interesting to do a Wilcott-Davis fit to see where the lanthanide is located here, but that takes us away from our main interest.

Dunning may also be defined as the act of making something dun colored-rather like this brief report.

Wishing you a better spring,

William B. Smith

56.4(0.08)25.9(0.14) **~131.5(0.35)** 150.0(0.22) CH3 0 29.7(0.16) 110.7(0.24) --40.9(0.22)108.4(0.46) 23.1(0.24) 43.4(0.35) 71.3(1.00)

55.6(0.73) 123.8(0.78)

CENTRO DE INVESTIGACIÓN DEL IPN

APARTADO POSTAL 14-740 MEXICO 14, D. F.

DEPARTAMENTO DE QUIMICA

April 20, 1983.

Professor Bernard L. Shapiro Department of Chemistry Texas A and M University College Station, Texas U. S. A. 77843

Fluorine-19 nmr measurements of steroidal drugs.

Dear Professor Shapiro:

The low solubility of many fluorinated corticosteroids in conventional nmr solvents has in general precluded systematic measurements of both their proton and fluorine spectra, although some isolated cases can be found in the literature (1).

Although the use of trichloroacetyl isocyanate reaction (2) of some corticosteroids allowed analyses of the derived carbamates, in particular for fluorine containing molecules a much more specific spectral determination is provided by F-19 nmr of the underivatized compounds.

Pulse-FT determinations could solve the sensitivity problem for these substances, although the very large spectral width required particularly for fluorine measurements preclude its application in a practical way.

In order to solve these problems spectra were measured using 12 mm sample tubes in which it was possible to dissolve some 5 mg of the individual steroids in approximately 4 ml of either deuterochloroform or hexadeuteroacetone. This concentration is still low for obtaining good fluorine spectra and therefore the measurements were performed by continuous wave time averaging using a Varian XL-100A spectrometer in combination with the 620L-100 computer. In general 10 to 25 scans were accumulated thus providing a 3 to 5 fold increase in sensitivity. The deuterium of the solvent was used to lock the spectrometer. The data show that there is a $\Delta\delta \simeq 20$ ppm between the 6α and the 9α fluorines in accordance with the known shielding of an axial as compared to an equatorial fluorine (3). Moreover, observation of the $^{19}F^{-1}H$ vicinal coupling constants evidences the potential of this technique for the conformational analysis of this type of molecules.

- 1.- A. D. Cross and P. W. Landis, J. Am. Chem. Soc., <u>86</u>, 4005 (1964).
- 2.- M. Lanouette, D. Legault and B. A. Lodge, <u>J. Pharm. Sci.</u>, <u>65</u>, 8, 1214 (1976).
- 3.- H. J. Schneider, W. Gschwendtner, D. Heiske, V. Hoppen and F. Thomas, Tetrahedron, 33, 1769 (1977).

I:
$$R_1 = F$$
; $R_2 = F$; $R_3 = H$.
II: $R_1 = F$; $R_2 = F$; $R_3 = Ac$.
III: $R_1 = F$; $R_2 = H$; $R_3 = H$.
IV: $R_1 = H$; $R_2 = F$; $R_3 = Ac$.

VII: R = H. VIII: R = Ac.

V:
$$R_1 = F$$
; $R_2 = F$; $R_3 = H$.
VI: $R_1 = H$; $R_2 = F$; $R_3 = Ac$.

IX

19 F CHEMICAL SHIFTS (in ppm)* AND HETERONUCLEAR 19 F-1 H COUPLING CONSTANTS (in Hz)** IN FLUORINATED CORTICOSTEROIDS

Compound	Solvent	<u>F</u> 6α				F _{9α}		
		δ	^J 6α,6β	^J 6α,7α	^J 6α,7β	δ	^J 9α,8β	^J 9α,11α
I	C ₂ D ₆ CO	186.8	48.5	14.5		165.5	28.0	9.0
	CDC1 ₃	188.1	49.0	13.5		166.7	27.0	8.3
II	C ₂ D ₆ CO	187.1	49.0	14.5		165.5	28.0	8.5
	CDC1 ₃	187.9	49.0	14.0		167.0	27.0	8.0
III	c ₂ D ₆ čo					165.9	29.0	10.0
	CDC1 ₃					167.6	28.0	9.0
IV	c ₂ D ₆ co	183.9	49.0	11.5	3.2			
	CDC1 ₃	185.0	49.0	11.7	3.5			
٧°	C ₂ D ₆ CO	186.8	48.5	14.5	-	164.2	28.0	9.5
VI	C ₂ D ₆ CO	185.8	50.0	13.5				
	CDC1 ₃	187.0	49.5	13.0				
VII°	C ₂ D ₆ CO		•			165.4	29.0	9.7
VIII	C ₂ D ₆ CO					165.4	29.7	10.0
	CDC1 ₃					167.0	29.0	9.0
IX	c ₂ D ₆ co	186.9	49.0	14.8		164.4	29.0	10.0
	CDC1 ₃	187.5	48.8	13.6		164.6	28.0	82
X	c ₂ D ₆ co	184.0	48.8	11.0	3.0			
	CDC1 ₃	184.9	49.0	12.0	3.5			•
XI	C ₂ D ₆ CO.	184.0	49.7	14.5	9.0		·	
	CDC1 ₃	185.2	49.7	14.5	9.5			

^{*}From CFC1 $_3$ as internal standard ($\delta \text{=0 ppm})\text{.}$

Sincerely yours,

P. Joseph-Nathan Professor of Chemistry R.L. Santillán

J. Espiñeira

^{**}Additional couplings are evident from signal broadening.

[°]Insoluble in CDC13 for measurement.

TEXAS A&M UNIVERSITY

DEPARTMENT OF CHEMISTRY COLLEGE STATION, TEXAS 77843

1983-84 Newsletter Subscriptions and Related Matters

It will soon be time for the execution of subscription renewals for the 1983-84 year of the TAMU NMR Newsletter - our 26th year of continuous operation! (I suppose that we will soon have to begin offering "Senior Citizen" discounts.) The mails being what they are, especially overseas postings, many of you will receive your subscription renewal notices and invoices about the same time, or even before, this May 1983 issue. In any event, permit me to make a special plea that subscription renewal handling be done expeditiously, for it really does add significantly to our cost of operation when we must do extra book-keeping, send reminders, etc. If you do not receive a subscription renewal notice and invoice by, say, July 18, please inform us immediately by airmail, and we will provide replacement copies. The subscription renewal notices will go into the mail on July 1 (airmail for overseas recipients).

I know it will not escape your attention that there is a small increase in the subscription fee for the Newsletter subscriptions, beginning with the 1983-84 year. The small increase only partially looks after the increased costs since the current rate went into effect in 1981.

The entirely self-supporting TAMU NMR Newsletter continues to exist, then, not solely by funds raised by the subscription fees we charge. Rather we continue to be dependent strongly on the generosity of our Contributors and Sponsors, and on the equally reliable and supportive role played by our Advertisers, most especially those who advertise every month. All of us in the NMR community owe these Advertisers a debt of gratitude for keeping this apparently useful Newsletter alive. Thus please permit me to ask that you provide suitable feedback to our Advertisers to the effect that the Newsletter is indeed worthy of their continued and generous financial support. We continue to need additional support in all categories, especially Sponsors, Contributors, and Advertisers.


I have in mind certain changes in the Newsletter - primarily additions - and will of course welcome suggestions for improving any aspect of the Newsletter operation. Please don't be shy! The first change which will occur will be the occasional, irregular appearance of a "Book Review" section. I have volunteered my good friend and colleague, Professor William B. Smith of Texas Christian University, to do the book reviews. Bill's breadth of interest and well-known fairness will make the book reviews useful; his incisiveness and wit guarantee that the reviews will also be a pleasure to read and occasionally amusing. (If the latter aspect begins to fail, we will apprentice him to Ray Freeman for a few months.)

Again, if you do not have a Newsletter subscription renewal notice and invoice in your hands by Monday, July 18, 1983, please let us know immediately. The Newsletter continues to thrive technically, but we must unfortunately keep its finances healthy as well...

Dany Rapico

B. L. Shapiro

2 May 1983



You never had NMR like this before.

NMR is the most powerful structural analysis technique available to the organic chemist today. It lets you determine molecular structure, investigate kinetic phenomena, and perform quantitative analysis of organic compounds in complex mixtures.

But most NMR systems are designed for research, not for quick and easy use.

That's why Nicolet developed the QE-300—a low cost, high performance NMR system for routine analysis of organic materials.

Optimized for analytical power

The QE-300 has all the power you need for high sample throughput, high production laboratory use. A 300 MHz superconducting magnet combined with optimized RF electronics and a specially designed single carbon/hydrogen probe provide you with greater chemical shift dispersion and sensitivity. So you can run more samples faster and interpret results with less ambiguity.

Complete computer control

The QE-300 makes it possible for virtually anyone to get high quality NMR spectra—quickly and easily—without, previous training or NMR experience.

A powerful NMR software package lets you set up and automate sample runs to your own criteria. And since the software is completely menudriven, it's almost impossible to make a mistake.

All an operator has to do is slip a sample in and type a single key on the control console. The computer controls the spin rate, magnet shimming, lock

frequency, and spectral phasing, acquires data to preset default parameters, and prints out the results, complete with full annotation of system settings.

You don't even have to change probes to get carbon and hydrogen spectra on the same sample. Just enter a single command, and the QE-300 makes the ¹H to ¹³C conversion automatically.

Attractive features The QE-300 has features that make it a pleasure to use, too. Like a color display, an eight-pen color plotter, and a dual floppy disc system for interactive spectral analysis and unlimited external data storage.

And just as importantly, the QE-300 won't lead you down an alley of obsolescence. Its open-ended design leaves plenty of room for adding even more high performance capability in the future.

The Nicolet QE-300. It's the ideal NMR system for routine analysis, QC, and troubleshooting.

Find out how friendly NMR can be. Circle the reader service number, or call Rich Bohn today at (415) 490-8300.

You never had NMR this easy.



255 Fourier Avenue, Fremont, CA 94539 (415) 490-8300 TWX 910-381-7025

THE GX SOLID EADER

With Multi-Nuclear/ **Multi-Field Solid State Probes**

High field solid sample probe for JEOL's 270 MHz SCM!

*Tunable heads — interchangeable plug-in matching units for observation of


13C (~67.8 MHz) 31P (~109.2 MHz) ²⁹Si (~53.6 MHz)

with one probe!

- Self-starting rotor/stator design!
- High Speed magic-angle double air bearing sample spinning (>4.0 KHz)!
- "Magic lift probe" for quick sample change and probe insertion!
- All this, in addition to a full line of dual and broad-band high resolution liquid sample probes!

SOLIDS UPDATE -NOW AVAILABLE

- MULTI-NUCLEAR SOLIDS PROBE FOR THE GX-270
- WIDE BORE MAGNET WITH VARIABLE TEMPERATURE SOLIDS FOR THE GX-270/89



235 Birchwood Ave., Cranford, NJ 07016

(201) 272-8820

TUNABLE

OBS COIL/STATOR

MAGNET MOUNTING BRACKET

PROBE BODY

TUNABLE