
Comments	. Shapiro, B. L.	2
40 Years in Show Business (<i>Second Installment</i>)	. Bothner-By, A. A.	5
Use of ^1H MAS NOESY to Determine Location and Diffusion in Small Molecules in Biomembranes	. Yau, W.-M., and Gawrisch, K.	7
Early Adventures of Ray Freeman , Master of NMR Choreography	. Anderson, W. A.	11
<i>A Blast from the Past</i> : Decomposition of Unresolved Multiplets	. Freeman, R., and Gestblom, B.	15
<i>Another Blast from the Past</i> : Early Jurassic Digital Computers	. Freeman, R.	17
<i>Still Another Blast from the Past</i> : Niveaux d'Energie A_nX_m	. Perchery, S.	18
HCN TROSY for Nucleic Acids	. Fiala, R., Czernek, J., and Sklenář, V.	21
Water Inversion Using a Radiation Damping Control Unit	. Bornet, G., Guerlesquin, F., Brevard, C., and Piotto, M.	23
Fast Lithium Sigmatropic Shift in an Otherwise Dynamically Frozen Allylic Lithium Compound	. Fraenkel, G., and Fleischer, R.	27
Equipment Available	. James, T. L.	28
Spectral Estimation of NMR Relaxation	. Naugler, D., and Cushley, R. J.	29
Spectral Estimation of NMR Relaxation	. Istratov, A.	30
Temperature Calibration of a Solid State NMR Probe	. Slupsky, C., and Sykes, B. D.	33
PSGE NMR Gradient Pulses and Artificial Diffraction Peaks	. Price, W. S.	35
Oxygen $T_{1\rho}$ Magnetic Relaxation Dispersion	. Teng, C.-L., Bryant, R. G., and Wagner, S.	37
Pseudorandom k -Space Sampling	. Bourgeois, M., Briguet, A., and Graveron-Demilly, D.	38
Glycine Revisited	. Taylor, R. E.	39

A monthly collection of informal private letters from laboratories involved with NMR spectroscopy. Information contained herein is solely for the use of the reader. Quotation of material from the Newsletter is *not* permitted, except by direct arrangement with the author of the letter, in which case the material quoted *must* be referred to as a "Private Communication". Results, findings, and opinions appearing in the Newsletter are solely the responsibility of the author(s). Reference to The NMR Newsletter or its previous names in the open literature is strictly forbidden.

These restrictions and policies apply equally to both the actual Newsletter recipient/participants and to all others who are allowed access to the Newsletter issues. Strict adherence to this policy is considered essential to the successful continuation of the Newsletter as an informal medium for the exchange of NMR-related information.

	Frequency Range	Resolution	Switching Time	Phase-Continuous Switching	Rack-Mount Cabinet Dim. ¹	Remote-Control Interface	Price Example ²
PTS 040	.1-40 MHz	optional .1 Hz to 100 KHz	1-20 μ s	optional	5¼"H×19"W	BCD (std) or GPIB (opt)	\$5,330.00 (1 Hz resol., OCXO freq. std.)
PTS 120	90-120 MHz	optional .1 Hz to 100 KHz	1-20 μ s	optional	5¼"H×19"W	BCD (std) or GPIB (opt)	\$5,330.00 (1 Hz resol., OCXO freq. std.)
PTS 160	.1-160 MHz	optional .1 Hz to 100 KHz	1-20 μ s	optional	5¼"H×19"W	BCD (std) or GPIB (opt)	\$6,495.00 (1 Hz resol., OCXO freq. std.)
PTS 250	1-250 MHz	optional .1 Hz to 100 KHz	1-20 μ s	optional	5¼"H×19"W	BCD (std) or GPIB (opt)	\$7,440.00 (1 Hz resol., OCXO freq. std.)
Type 1 PTS 310 Type 2	.1-310 MHz	1 Hz	1-20 μ s	standard	3½"H×19"W	BCD (std) or GPIB (opt)	1 Hz resol., OCXO: \$6,425.00 1 Hz resol., OCXO: \$5,850.00
PTS 500	1-500 MHz	optional .1 Hz to 100 KHz	1-20 μ s	optional	5¼"H×19"W	BCD (std) or GPIB (opt)	\$8,720.00 (1 Hz resol., OCXO freq. std.)
PTS 620	1-620 MHz	optional .1 Hz to 100 KHz	1-20 μ s	optional	5¼"H×19"W	BCD (std) or GPIB (opt)	\$9,625.00 (1 Hz resol., OCXO freq. std.)
PTS 1000	0.1-1000 MHz	optional .1 Hz to 100 KHz	5-10 μ s	optional	5¼"H×19"W	BCD (std) or GPIB (opt)	\$11,830.00 (1 Hz resol., OCXO freq. std.)
PTS 3200	1-3200 MHz	1 Hz	1-20 μ s	optional	5¼"H×19"W	BCD (std) or GPIB (opt)	\$14,850.00 (1 Hz resol., OCXO freq. std.)
PTS x10	user specified 10 MHz decade	1 Hz	1-5 μ s	standard	3½"H×19"W	BCD (std) or GPIB (opt)	\$3,000.00 (1 Hz resol., OCXO freq. std.)
PTS D310	two channels .1-310 MHz	.1 Hz	1-20 μ s	standard	5¼"H×19"W	BCD (std) or GPIB (opt)	\$8,560.00 (1 Hz resol., OCXO freq. std.)
PTS D620	two channels 1-620 MHz	.1 Hz/.2 Hz	1-20 μ s	standard	5¼"H×19"W	BCD (std) or GPIB (opt)	\$13,240.00 (.1 Hz/.2 Hz resol., OCXO freq. std.)



1 Bench cabinets are 17" wide.

2 Prices are U.S. only and include Manual and Remote (BCD) Control; PTS 3200 Digital Front Panel.

PTS CAN SUPPLY OEM-TYPE SYNTHESIZERS FOR ALL LEADING NMR-SPECTROMETER PRODUCTS.

PROGRAMMED TEST SOURCES, INC.

P.O. Box 517, 9 Beaver Brook Rd., Littleton, MA 01460 Tel: 978-486-3400 Fax: 978-486-4495

http://www.programmedtest.com • e-mail: sales@programmedtest.com

THE NMR NEWSLETTER		NO. 489, JUNE 1999		AUTHOR INDEX	
Anderson, W. A.	11	Czernek, J.	21	Guerlesquin, F.	23
Bornet, G.	23	Fiala, R.	21	Istratov, A.	30
Bothner-By, A. A.	5	Fleischer, R.	27	James, T. L.	28
Bourgeois, M.	38	Fraenkel, G.	27	Naugler, D.	29
Brevard, C.	23	Freeman, R.	15, 17	Perchery, S.	18
Briguet, A.	38	Gawrisch, K.	7	Piotto, M.	23
Bryant, R. G.	37	Gestblom, B.	15	Price, W. S.	35
Cushley, R. J.	29	Grav'n.-Dem'y. D.	38	Shapiro, B. L.	2
				Sklenář, V.	21
				Slupsky, C.	33
				Sykes, B. D.	33
				Taylor, R. E.	39
				Teng, C.-L.	37
				Wagner, S.	37
				Yau, W.-M.	7

THE NMR NEWSLETTER		NO. 489, JUNE 1999		ADVERTISER INDEX	
Advanced Chemistry Development, Inc.	25	JEOL	outside back cover		
Aldrich Chemical Company, Inc.	31	Programmed Test Sources, Inc.	inside front cover		
AMT	9	Varian, Inc.	3		
Bruker Instruments, Inc.	19				

SPONSORS OF THE NMR NEWSLETTER

Abbott Laboratories	JEOL (U.S.A.) Inc., Analytical Instruments Division
Aldrich Chemical Company, Inc.	The Lilly Research Laboratories, Eli Lilly & Company
AMT	Merck Research Laboratories
Amgen, Inc.	Nalorac Corporation
Anasazi Instruments, Inc.	Oxford Instruments
Astra AB	Pharmacia & Upjohn, Inc.
Avanti Polar Lipids, Inc.	Programmed Test Sources, Inc.
Bruker Instruments, Inc.	SINTEF Unimed MR Center, Trondheim, Norway
Cambridge Isotope Laboratories	Tecmag
Cryomag Services, Inc.	Unilever Research
The Dow Chemical Company	Union Carbide Corporation
E. I. du Pont de Nemours & Company	Varian, Inc.
Hewlett-Packard Company	Zeneca Inc.
Isotec, Inc.	

FORTHCOMING NMR MEETINGS

Gordon Conference on Magnetic Resonance, Henniker, NH, **June 27 - July 2, 1999**. Contact: Robert Tycko, tycko@helix.nih.gov. See Newsletter 487, 37.

Royal Society of Chemistry: 14th International Meeting on NMR Spectroscopy, Edinburgh, Scotland, **June 27 - July 3, 1999**; Contact: '99NMR14' c/o Mrs. Paula Whelan, The Royal Society of Chemistry, Burlington House, London W1V 0BN, England; +44 0171 440 3316; Email: conferences@rsc.org\

Rocky Mountain Conference NMR Symposium, Denver, CO, **August 1-5, 1999**. See Newsletter 487, 36 and <http://india.cchem.berkeley.edu/~rmc/>.

SMASH No. 1 (Small Molecules Are Still Hot), Argonne, IL, **August 15-18, 1999**; Contact: Ms. Karen McCune, (mccune_karen_a@lilly.com, 317-276-9783) or S. R. Maple (maple_steven_r@lilly.com) or G.E.Martin (gary.e.martin@am.pnu.com) or A. G. Swanson (alistair_swanson@sandwich.pfizer.com). See Newsletter 487, 17.

Applications of NMR to Complex Systems, symposium at the American Chemical Society Meeting, New Orleans, LA, **August 22-26, 1999**; Contact: R. E. Botto, Symposium Chair, Chemistry Division, Argonne National Laboratory, 9700 South Cass Ave., Argonne, IL 60439; 630-252-3524; Fax: 630-252-9288; E-mail: robert_botto@qmgate.anl.gov

Medical Imaging: NMR and Nuclear Tracers, colloquium at the 12th Entretiens Jacques Cartier, Lyon, France, **December 5-8, 1999**; See <http://jade.univ-lyon1.fr/JacquesCartier/> and Newsletter 488, 38.

41st ENC (Experimental NMR Conference), Asilomar Conference Center, Pacific Grove, CA, **April 9-14, 2000**; Contact: ENC, 1201 Don Diego Avenue, Santa Fe, NM 87505; (505) 989-4573; Fax: (505) 989-1073; E-mail: enc@enc-conference.org.

XIX International Conference on Mag. Res. in Biological Systems, Florence, Italy, **August 20-25, 2000**. Contact: Profs. Ivano Bentini or Lucia Banci, Chem. Dept., Univ. of Florence, Via G. Capponi 7, I-50121, Florence, Italy; Phone: +39-055-2757600; Email: icmrbs@lrm.fi.cnr.it; Fax: +39-055-2757555; <http://www.lrm.fi.cnr.it/icmrbs.html>.

Additional listings of meetings, etc., are invited.

THE
NMR
NEWSLETTER

B. L. Shapiro, Publisher
 The NMR Newsletter
 966 Elsinore Court
 Palo Alto, CA 94303-3410

650-493-5971
 650-493-1348 Fax
 shapiro@nmrnewsletter.com
 http://www.nmrnewsletter.com

Comments.

1. *"Spin Choreography, a Symposium in Appreciation of Ray Freeman"*. Starting on p. 11 with a narrative by **Wes Anderson** are some pages relevant to the April symposium at the University of Cambridge. The Newsletter also has reasons to appreciate Ray, for he holds the record - by a wide margin - for the number of technical contributions by one person in the Newsletter's pages: from his first missive in June of 1961 (issue no 33) to the present, Ray has authored no fewer than 68 letters! Several of these have featured not only first-rate science, but also writing distinguished by its literacy and wit, and also by screamingly funny and appropriate cartoons. For all of these aspects, and most especially for his loyalty over so many years, we are most grateful.

If any other participants (or indeed non-participants) in the Cambridge symposium would like to share thoughts/experiences/photos/cartoons/whatever concerning Ray, we would be happy to receive such pages (camera-ready, svp) for the July or August Newsletter issues (deadline dates are June 25 and July 23).

2. *"40 Years in Show Business"*, by **Aksel Bothner-By**. This imaginatively formatted saga continues on p. 5. Although I am grateful for these efforts, allow me to correct the impression that the inception of the Newsletter in October of 1958 was solely my doing - my clear recollection was that it was a joint effort until I left Mellon Institute early in 1964. The first name of the newsletter, MELLONMR (Monthly Ecumenical Letters from Laboratories of NMR, pronounced 'mellonmer') was Aksel's doing. The operating ideas he and I settled on for the Newsletter have continued to this day.


3. *"Blast from the Past"*. Are these of interest to anyone other than the authors and the publisher? Any suggestions for future examples.

4. Some small practical matters:

- a) Another plea to include in your technical letters your phone and fax numbers, and especially, your **e-mail addresses**, to facilitate communication.
- b) Subscription invoices and requests for Sponsorship renewals will be mailed on July 1, for your kind attention. Subscription rates will again remain unchanged.
- c) Reader feedback is always appreciated - even when I am unable to agree or acquiesce!

Barry Shapiro
 June 1, 1999.

Experience Experimental Breakthroughs




With Varian's gHX Nano Probe

The pioneer and leader in high resolution magic angle spinning technology, Varian continues to break sensitivity barriers with its new gHX Nano™ Probe. For the first time, you can achieve unmatched experimental performance by combining indirect detection capability with "Magic MAS" gradients and the unique Nano Probe technology. The gHX Nano Probe is the only probe with "Magic MAS" gradients, highest sensitivity proton detection, and a broadband X nucleus decoupling coil, making it a powerful tool for even the most difficult sample.

The probe is ideal for on-bead chemical applications and solvent suppression experiments where dual channels and indirect detection are necessary.

See for yourself how Varian can help you with your combinatorial, metabolite, or natural product applications. For more information on any of our NMR probes or accessories, call your local sales office or check out the web at www.nmr.varian.com.




- 1992 Varian's Nano Probe changes the way natural product chemists analyze small samples by breaking the sensitivity barrier for high resolution NMR.
- 1993 Varian provides undreamed of sensitivities for obtaining ¹³C NMR spectra for tiny amounts of solution-state samples.
- 1994 Varian brings high resolution NMR to chemists analyzing compounds still bound to solid-phase synthesis resins with the second generation Nano Probe. This capability has since blossomed into the rapidly growing field of HR-MAS, whereby a wide range of heterogeneous samples are yielding to NMR analysis.
- 1996 Varian publishes the definitive "High-Resolution NMR Spectra of Solid-Phase Synthesis Resins."
- 1997 Varian announces the innovative NMS, bringing automation to Nano Probes.
- 1998 Varian announces the gHX Nano Probe which goes even further by offering indirect detection capabilities for even the most demanding experiments.

Manufacturing Facilities Varian NMR Instruments, Building 4, 3120 Hansen Way, Palo Alto, California 94304-1030, Tel 650.424.4876, Fax 650.852.9688
 • <http://www.nmr.varian.com> • **Argentina** Buenos Aires, Tel 1.783.5306 • **Australia** Mulgrave, Victoria, Tel 3.9566.1133 • **Brazil** São Paulo, Tel 11.820.0444
 • **Canada** Mississauga, Ontario, Tel 1.800.387.2216 • **France** Les Ulis, Tel 1.69.86.38.38 • **Germany** Darmstadt, Tel 06151.7030 • **India** Mumbai, Tel 22.859.0178
 • **Italy** Milan, Tel 2.921.351 • **Japan** Tokyo, Tel 3.5232.1211 • **Korea** Seoul, Tel 2.3452.2452 • **Mexico** Mexico City, Tel 5.523.9465 • **Netherlands** Houten, Tel 3063.50909 • **Switzerland** Basel, Tel 61.295.8000 • **Taiwan** Taipei Hsien, Tel 2.698.9555 • **United Kingdom** Walton-on-Thames, Tel 1932.898.000
 • **United States** California, Tel 800.356.4437 • **Venezuela** Valencia, Tel 4125.7608
 • Other sales offices and dealers throughout the world

40 YEARS IN SHOW BUSINESS

Dear Barry- The adventures continued (2)

We wrote to G. Traverso, in Italy, and told him we could decide between two structures for him



A year later , X-ray diffraction




S. Bezzi, M. Mammi, & C. Garbuglio

Holy Shinola!


Aksel moves to Mellon Institute, in Pittsburgh 1958



Barry Shapiro and Cecile Maar came to Pittsburgh, and helped set up the new 60 MHz spectrometer.



Cecile shims the 60. For the meaning of "chocolat", consult a Belgian friend.



CHOCOLAT!

1958 began another 40 year show:
In October Barry brought out

Monthly
Ecumenical
Letters from
Laboratories
Of
N-M-R


No.1

(received 5/20/99)

40 YEARS IN SHOW BUSINESS

Cleveland

1960 - Bill Ritchey organizes the first ENC, another 40 year travelling show.




Barry attends., and his enthusiasm convinces everybody we should hold the next ENC at Mellon Institute. Turned out to be ten ENCs



SMASH HIT

2nd ENC

A60
the
new
INTER




With a lot of head-scratching and poring over John Pople's chapters in P.S.&B., we wrote FREQINT, to compute line freqs and intensities from J 's and δ 's

3rd ENC



**We got a big
boost in 1962
when Turi
Castellano joined
us.**



5th ENC -

LAOCOON - to get J 's and δ 's by fitting the line freqs. Lines often crossed over when J or δ was varied, creating snakes.



Andrew C. Boothby
ab62@andrew.cmu.edu



Dr. B.L. Shapiro
The NMR Newsletter
966 Elsinore Court
Palo Alto, CA 94303

National Institutes of Health
Bethesda, Maryland 20892

(received 5/13/99)
May 10, 1999

Use of ^1H MAS NOESY to determine location and diffusion of small molecules in biomembranes

Dear Dr. Shapiro:

The location of small molecules in biomembranes, when combined with data on partitioning and interaction enthalpies, provides very valuable information about the nature of interaction of these molecules with biomembranes. We have investigated location of Trp analogs [1], alcohol [2] and benzene [in preparation] by two-dimensional nuclear Overhauser enhancement spectroscopy (NOESY) on biomembranes. High resolution conditions have been achieved by magic angle spinning (MAS).

For the quantitative analysis of cross-relaxation rates, the NOESY diagonal peak and crosspeak intensities are recorded as a function of mixing time. At short mixing times, each crosspeak and the corresponding diagonal peak of the lipid resonance are treated as a pair of interacting protons [3]. The corresponding cross-relaxation rate, σ_{ij} , has been calculated by fitting crosspeak volumes to the equation

$$A_{ij}(t_m) = (A_{ij}(0)/2)(1 - \exp(-2\sigma_{ij}t_m))\exp(-t_m/T_{ij})$$


where $A_{ij}(t_m)$ is the respective crosspeak volume at mixing time t_m , $A_{ij}(0)$ is the diagonal peak volume at mixing time zero, and $1/T_{ij}$ is the rate of magnetization leakage towards the lattice.

For both aromatic compounds and alcohol, the resulting per-proton cross-relaxation rates reflect a preferential interaction of these molecules with the membrane/water interface region. At the same time, there is weaker cross-relaxation to resonances that are located in the bilayer hydrophobic center, like the terminal methyl protons of hydrocarbon chains. We conclude that location of such molecules in the bilayer is more accurately described by a distribution function along the bilayer normal. Furthermore, the wide range of contacts is also a reflection of the tremendous disorder in lipid packing of a liquid-crystalline lipid matrix [2-4]. For all investigated substances, we were able to verify the preference for a biomembrane interface location by other spectroscopic means. We have strong evidence that the averaged magnitude of cross-relaxation rates depends on the diffusion rates of the small molecules in the lipid matrix. This allows new insights into diffusion of small molecules in the plane of the lipid matrix and across the bilayer. For example, the tryptophan analogs appear to have diffusion rates that are comparable to the lateral lipid

diffusion rates, while indene, benzene, and alcohol diffuse much more rapidly within the bilayer.

- [1] Wai-Ming Yau, William C. Wimley, Klaus Gawrisch, and Stephen H. White, (1998) *Biochemistry* **37**, 14713-14718.
 [2] Laura L. Holte and Klaus Gawrisch, (1997) *Biochemistry* **36**, 4669-4674.
 [3] Daniel Huster, Klaus Arnold, and Klaus Gawrisch (1999) *Journal of Physical Chemistry* **103**, 243-251.
 [4] Daniel Huster, Klaus Gawrisch (1999) *J. Am. Chem. Soc.* **121**, 1992-1993.

Sincerely yours,



Wai-Ming Yau




Klaus Gawrisch

Section of NMR
 Laboratory of Membrane Biochemistry and Biophysics
 National Institute on Alcohol Abuse and Alcoholism

Phone: (301) 594-3750

Fax: (301) 594-0035

e-mail: gawrisch@helix.nih.gov



Per-proton normalized cross-relaxation rates between selected lipid resonances and membrane-incorporated aromatic molecules.




SCIENTIFIC & MEDICAL PRODUCTS



Model 3205 - 6 MHz to 220 MHz, 300 W, NMR Amplifier



Model 3445 - 10 MHz to 130 MHz, 2.0 kW, MRI Amplifier



200 MHz to 500 MHz, 50 W, NMR Module



6 MHz to 220 MHz, 300 W, NMR Module



*Model 4T70 - 25 MHz to 175 MHz,
7.0 kW, MRI Amplifier*

AMT's scientific products are used extensively in Nuclear Magnetic Resonance (NMR) systems. These amplifiers cover the frequency ranges of 6 MHz to 950 MHz, with power levels as high as 2.0 kW peak power at 10% duty cycle.

AMT's medical products are employed in Magnetic Resonance Imaging (MRI) systems. These amplifiers cover the frequency ranges of 10 MHz to 200 MHz with power levels as high as 12.0 kW peak power at 10% duty cycle.

All amplifiers have dual mode capability and can be operated in either a pulsed or CW mode. Scientific and Medical customers include both OEM system manufacturers and end users.



COMPANY

AMT designs, develops and manufactures custom radio frequency (RF) and microwave power amplifiers for the wireless, scientific/medical and application specific industries. The company has been in business since 1984 and currently has over 60 employees, including 20 experienced engineers.

AMT has a worldwide reputation as a leading supplier of high power, solid state power amplifier products that operate at frequencies between 1 MHz and 3 GHz and provide RF power from several watts to several kilowatts. Its products are noted for their exceptional performance, highest quality and superior reliability.

The company's products are sold to numerous major corporations, universities and research centers throughout the world.

FACILITIES

AMT is located in Anaheim, California and occupies a 25,000 square foot facility allocated to engineering, manufacturing, quality assurance, marketing/sales, administration and finance.

Engineering areas include an R & D laboratory, a tool and die shop, mechanical design and drafting areas, an environmental testing laboratory and document control. The R & D laboratory is equipped with all of the latest design and testing equipment including intermodulation distortion simulators, network analyzers, spectrum analyzers, signal generators, noise figure meters and infrared (IR) scanners. The environmental testing laboratory includes equipment to simulate shock, vibration and thermal environments.

Manufacturing areas include a controlled access stock room, a 10,000 square foot assembly area and a production test area employing automatic testing. Also included is an environmental laboratory used for environmental stress screening of production products.

PRODUCTS

AMT's products vary in complexity from single modules, to rack-mounted amplifiers, to complete transmitter systems. The rack-mounted amplifiers and complete transmitter systems typically include detection/protection circuitry, built-in power supplies, front panel metering and digital and/or analog interface controls. Both forced air and/or water cooling are used, depending on the customer's requirements.

AMT's products feature highly reliable technical solutions designed for producibility and reliability. Producibility is enhanced through the use of surface mount components and circuit designs that eliminate the need for excessive alignment during the production cycle. High reliability is accomplished through the implementation of conservative thermal and RF circuit design and sophisticated self-protection schemes. Reliability is further enhanced during the design phase by employing detailed environmental testing.

These factors, along with computer driven automatic testing and environmental stress screening of the final product, ensure that the performance, quality and reliability meet AMT's exacting standards.



An Employee Owned Company

2570 East Cerritos Avenue, Anaheim, CA 92806 • Tel: (714) 456-0777 • Call Jerry Toll Free: (888) 545-4AMT • Fax: (714) 456-0778

Early Adventures of Ray Freeman, Master of NMR Choreography

(as seen from Palo Alto)*

by

Weston A. Anderson

Varian, Inc.

This narrative* relates a selection of the many contributions Ray has made to NMR, and includes some of the interesting and sometimes humorous events that happened early in his career. It traces his path from his high school graduation up through his employment at Varian.

Ray was born in Derbyshire on January 6, 1932, and attended Nottingham High School, graduating in 1949. From high school, Ray went directly into the Air Force for a year and a half, and then went up to Oxford University to study chemistry, receiving his Bachelor of Science degree in 1955.

At Oxford he got his first taste of NMR, working in Rex Richards' group. As one of his projects, he prepared a sample of potassium amide and measured the spectrum of the solid. (He used the calculations of Pake to determine the inter-proton distance in the amide ion). This may have been more challenging than it appears, as Richards has recently commented, "Making the stuff without blowing yourself up is quite difficult, and nowadays would only be allowed with extreme precautions!" Ray's response: "All the compounds Rex suggested were dangerously explosive, or poisonous, or both. After the NMR measurements, we buried the samples in a grassy slope just outside the lab." As part of his D. Phil. work, Ray built a Pound spectrometer to measure the chemical shifts of compounds of gallium, thallium and cobalt, his first experience measuring chemical shifts in solution. He received his D. Phil. Diploma in 1957.

Saclay

After his D. Phil. Degree from Oxford, Ray went to Saclay, France, to work in Anatole Abragam's group. There he met Ionel Solomon and Bob Pound, who was on sabbatical from Harvard. At that time Pound was concerned with radiation damping. (The marginal oscillators used by Pound and others enhanced the radiation damping, causing line broadening of all high-resolution NMR lines). Based on some of Pound's ideas, Ray went to work to build a high-resolution spectrometer which included an internal lock to maintain field/frequency stability, and audio frequency field modulation and sideband detection to provide base line stability. This super-regenerative oscillator had the nice property that it would automatically lock onto an NMR line.

During Ray's time at Saclay, Prof. Abragam is known to have quipped: "In the midst of the broken English that is heard here, Ray, you remain our standard of English" to which Ray replied: "If I am to serve as a standard, it is high time I was sent home to be recalibrated". In French, the word for "standard" is *étalon*, which also means "stallion", so Ray suspected Abragam had a pun in mind!

*Excerpted from the talk presented by W.A.A. on April, 8, 1999 at the University of Cambridge in the "Spin Choreography, a Symposium in Appreciation of Ray Freeman."

Ray made several trips to Paris, during one of which he remembers a meeting where Varian's Jim Shoolery and Emery Rogers gave a talk announcing the Varian HR-60 NMR spectrometer. This was a great shock to Ray, who in 1957 was still working at 30 MHz.

National Physical Laboratory

In 1959, Ray went to the National Physical Laboratory at Teddington, England, where he became enchanted with NMR double resonance. He added a sideband double resonance device to NPL's HR-60 spectrometer. Later he added an additional sideband channel to provide a lock channel. Dennis Evans, a colleague from his graduate student days at Oxford, provided an interesting application. He had studied the proton-thallium spin couplings in the diethylthallium cation. There were two methyl groups and two methylene groups, both coupled to the thallium spin. Since the effects of selective irradiation of the thallium frequencies depended upon the magnitude and signs of the two thallium-proton coupling constants, Evans was able to determine the relative signs of the thallium-proton couplings. When Evans described this experiment, Freeman and David Whiffen realized that similar measurements could be carried out on any three-spin system, provided one could obtain adequate sensitivity. They demonstrated that in 2,3-dibromopropionic acid*, the *gem*- and *vic*- couplings have opposite signs.

Varian Sabbatical

In 1961, Ray took a sabbatical leave to work at Varian. He arrived in November, with his wife, Anne Marie, and two children, Dominique and Anne. One of the things that Ray had to contend with was the fact that there is "the English language" and there is "American English". Ray soon found out that 'Berkeley' is not pronounced 'Barclay'. It seems there is a town in Somerset called Berkeley but pronounced Barclay. There were a number of other words that apparently were not in Ray's vocabulary. Some of these became apparent when Ray purchased a house from another Varian employee. When inspecting the house Ray found the Telly was called a TV, the Tap was called a Faucet, the Gar'-rage was called a Garage'. Being a little bit intimidated Ray pointed to a hoe and asked what that was called, and found it was a hoe, and his friend asked him "Don't you have them in England?"

There were also the American holidays. One day Ray arrived at work at the usual time and found the parking lot empty and the doors locked. That was the time that he learned that Washington's Birthday was a holiday.

One of the activities of the R&D personnel in those days was an annual back-packing trip to the mountains. We would choose a scenic spot in the Sierra Nevada mountain range and then would backpack in for 4 or 5 days. We would carry all of our own food and supplies, and usually little tents made of plastic, as it often rains in late afternoon and at night. Ray heard about our plans and was eager to go. Although it was a hard climb, there were many scenic views along the way. However the next morning we found that one person had a hard time crawling out of his sleeping bag. I still remember Ray's words: "Leave me alone!"

At the end of Ray's sabbatical year he was offered a permanent job at Varian. He accepted, taking an eight-month leave of absence from Varian before starting.

*It is well known that Ray holds the world's record for publishing papers using this cheap, readily available compound.

Ray's Scientific Contributions while at Varian

While at Varian Ray made many contributions to the NMR field. His main emphasis was the development of techniques for analyzing complex NMR spectra. One of the questions that came out of some of our informal seminars was the idea we called "spin tickling": What happens when a second RF* field affects just one transition, *i.e.*, if that RF field is too weak to cause decoupling?

Ray and I analyzed the basic idea, and showed that if one applies a weak perturbing RF frequency centered on one line, then all lines that had an energy level in common with the irradiated line would be split into doublets. This enabled one to trace out complex spin energy level diagrams. Furthermore, one could distinguish whether the two transitions were between levels with the same total magnetic quantum number M , or whether the quantum number of one of the levels differed by 2 units. These were called 'regressive' or 'progressive' levels, respectively.

Ferretti and Freeman illustrated a variation of this technique. They centered an observing RF on one line and then swept a second RF field through the entire spectrum. An amplitude change followed by a Torrey oscillation was observed when the second RF field passed through a transition line that had a common energy level. In addition, one could distinguish between progressive or regressive arrangements by whether there was an initial increase or decrease of the observed line intensity.

Another interesting spin decoupling technique we called 'modulation through the spin coupling'. Ray was interested in exploring the technique to see if it could be useful to solve real problems. We showed how it could be useful to detect carbon-13 satellites in proton spectra that might be hidden under the much stronger signals of protons coupled to carbon-12.

Ray was also quick to pick up on the technique of Fourier Transform NMR. An important question in the early days of FT NMR was whether one could achieve a high resolution NMR spectrum with the limited memory sizes that were available on computers at that time. Freeman and Jones showed they could pick up some very small long-range couplings of the aldehyde resonances in 3-bromothiophene-2-aldehyde. Here they used the concept of 'zero filling' the data, to improve the definition with a limited number of data points.

Another question was if relaxation times could be measured. Freeman and Hill combined the well-known techniques developed by Hahn and others into a series of Fourier Transform experiments, to measure all of the relaxation times T_1 and T_2 in a spectrum more or less simultaneously. For T_1 measurements, a series of 180° - 90° pulses were used, and a spin-lock sequence was used to measure T_2 .

*RF = radio frequency or Ray Freeman or République Française or *rinforzando*? (Ed.)


Oxford and Cambridge

In October 1973, Ray returned to England to accept a lectureship at Oxford, and in 1987, he moved to Cambridge as professor of chemistry. Here he and his coworkers developed many innovative pulse sequences, easily recognized by their descriptive and clever acronyms including BURP, ψ -COSY, CUSPIDOR, DANTE, EXORCYCLE, GARP, GROPE, INADEQUATE, INEPT, MLEV, SLURP, SPITTOON, TANGO, TSETSE, WALTZ, WURST, etc.

Ray has published 2 books and over 200 scientific papers. Last July, the University of Durham gave him an honorary degree.

Today, then, we salute **RAY FREEMAN, MASTER OF SPIN CHOREOGRAPHY.**

Palo, Alto, California
April 1999.



(1998) Self-captioned. Those who were not at the Cambridge symposium may try to identify Ray's companion - the first correct answer gets a prize.



(1958) Starting in September 1999, bookings for this stellar musical attraction are available for wakes, weddings, bar mitzvahs, etc. Apply to pescadero@marsals.com

A Blast from the Past

107-60

17 July 1967



varian/611 hansen way/palo alto/california 94303/u.s.a./415/326-4000

Professor B. L. Shapiro
Department of Chemistry
Stanford University
California, 94305


Decomposition of Unresolved Multiplets

Dear Barry,

For some time we have been interested in the problem of generating "differential" double resonance spectra, where all resonances are cancelled except for those perturbed by the irradiation. Experimentally this entails pulse modulation of H_2 at a very low frequency, and inversion of the signal in synchronism with the pulses, either in a conventional lock-in detector or through the "add" and "subtract" functions of a time averaging computer. Although this proved to be straightforward for heteronuclear double resonance¹, our first attempts at proton-proton work could only be described in terms of a phenomenon first reported by Thurber's grandmother², i.e. "modulation was leaking out all over the molecule." Every transition appeared to be directly connected to every other. We feel this was due to the rapid transients induced in the internal lock field-frequency loop as H_2 was switched, displacing all the spectral lines. It can be cured by introducing a short blanking period during which both the "add" and "subtract" functions are inhibited:




Among the several possible applications we would like to describe the decomposition of an unresolved multiplet--a 0.05 Hz doublet in the spectrum of 3-bromothiophene-2-aldehyde due to long-range coupling between the proton at position 4 and the aldehyde proton. This doublet consists of transitions X1 and X2, distinguishable by the different ways in which they fit into the energy level diagram; X1 is connected to M1 while X2 is connected to M3. Thus weak irradiation of line M1 causes X1 to split into a doublet but leaves X2 unaffected. Pulse modulation of the irradiating field on M1, followed by synchronous detection of the signal results in a response with a central



positive peak (X1) flanked by two weaker, negative signals (X1 split and inverted) with negligible contributions from X2.

Our time averaging device was a Mnemotron CAT 400 which has the feature that signals in alternate channels may be steered into two quite separate regions of the memory (let us call them the odd and even locations). At read-out, first all the odd channels are displayed sequentially and then, starting again at the left-hand edge of the recorder, the even channels. We therefore alternately irradiated lines M1 and M3, storing the X1 response in the odd locations and the X2 response in the even locations. The resulting trace is shown in the Figure.



Transient nutations of the magnetization vectors persist for a considerable time after H_2 is switched. It should be pointed out that there are certain features about these transient effects that are incompletely understood--in particular, oscillatory signals could be detected at a frequency that appeared to depend on the levels of both H_1 and H_2 . Unless the pulse period was increased sufficiently (22 secs) to minimize such transient effects, there was a detectable amount of "cross-modulation" between the X1 and X2 responses.

Yours sincerely,

Ray

Bo Gestblom

Ray Freeman

Bo Gestblom³

¹J. Chem. Phys. 43, 3087 (1965).

²J. Thurber, My Life and Hard Times, Chapter II.

³Varian Postdoctoral Fellow, on leave from the Institute of Physics, Uppsala.




Another Blast from the Past

UNIVERSITY OF CAMBRIDGE
DEPARTMENT OF CHEMISTRY
Lensfield Road, Cambridge, U.K.

Professor Barry Shapiro,
Editor, TAMU Newsletter,
Elsinor Court,
Palo Alto,
California 94303 USA

6 . ix . 93
(received XI/C9/XCIII)



Dear Barry,

Early Jurassic Digital Computers

I never quite understood how one manages to do long division on an abacus but I can imagine it can cause real difficulties when the denominator is zero. There were similar problems with digital computers in the 1960's; they seemed to blow all their fuses in this situation. Later computers checked one's arithmetic beforehand and sent a rather shirty error message. Division by zero was a crime almost as unforgivable as using double spacing in the Newsletter. *

Now there is a nice experiment suggested by Axel Bothner-By and Joe Dadok¹ where they deconvolute J-splittings from individual spin multiplets by dividing the corresponding time-domain signal by $\cos(\pi J^*t)$ where J^* is a trial coupling constant. Problems arise when the cosine wave is zero or near zero. Because of these numerical instabilities, a considerable amount of research has gone into alternative methods of deconvolution that avoid this problem. For example, we once advocated a method that carefully sampled the cosine wave in such a manner as to avoid the zero-crossings.²

We have recently attacked this problem again from another angle. In this particular case the cosine wave is a noisy signal which makes the problem even trickier. We test each ordinate of the time-domain denominator to see if its modulus is less than a predefined threshold θ . If it is, it is replaced by θ , retaining the original sign. A good value of θ seems to be something comparable with the peak-to-peak noise level, but in fact the method is surprisingly tolerant of the choice of θ (over several orders of magnitude). It even works when θ is so high that all the ordinates are replaced by $+\theta$ or $-\theta$ although in this case there are complications due to harmonics of the basic frequency. After division, any residual spikes in the quotient can be smoothed by a simple 1:2:1 convolution algorithm, repeated if necessary. We think that this might be one of the simplest solutions to the J-deconvolution problem.

Kindest regards,

Roy

Ray Freeman

¹ A. A. Bothner-By and J. Dadok, *J. Magn. Reson.* 72, 540 (1987).

² J. M. Le Parco, L. McIntyre and R. Freeman, *J. Magn. Reson.* 97, 553 (1992).

* Or misspelling Aksel, or leaving the last letter off Elsinore.

Thanks, Roy.

BLS.

Still Another Blast from the Past

100- 24

Académie des Sciences de Tahiti
Service de la Physique Expérimentale
et Résonance Magnétique
Avenue Louis de Bougainville, 13
Papeete

le 4 janvier 1967

M. le Professeur B. L. SHAPLIERO
Department of Chemistry
Illinois Institute of Technology
Chicago, Illinois 60616 U.S.A.

Mon cher Barrie,

Lors de votre récent passage à notre laboratoire, vous avez sollicité une communication digne de la centième édition de l'I.I.T.N.M.R.

Depuis quelque temps nous étudions les signes relatifs des constantes d'interaction spin-spin à longue distance dans certains extraits de l'huile de coprah. Au cours de ce travail nous avons éprouvé quelques difficultés à établir le réseau des niveaux d'énergie pour certains systèmes de deux groupes de spins nucléaires, par exemple A_2X_3 (selon la nomenclature de POPLS). Le problème a été résolu par une nouvelle méthode de construction graphique qui s'explique par l'exemple suivant.

On trace d'abord sur le graphique (Fig. 1) le profil du spectre A_2X_3 , se rappelant qu'un quadruplet 1:3:3:1 doit être décomposé en quadruplet (1:1:1:1) auquel on superpose un doublet (2:2) par dessus les deux lignes centrales, tandis qu'un triplet (1:2:1) est décomposé en triplet (1:1:1) avec un singulet central. Ainsi chaque résonance est en réalité deux sous-spectres¹ qui doivent être traités séparément.

Avec un compas centré sur chaque ligne d'un sous-spectre A quelconque, on trace des arcs qui se coupent avec d'autres arcs (à rayon r égal) centrés sur les lignes d'un sous-spectre X. Les intersections de ces arcs définissent les niveaux d'énergie, et, dans ce cas particulier de deux sous-spectres A et deux sous-spectres X, établissent quatre systèmes indépendants de niveaux (voir la figure ci-contre). Les transitions entre niveaux d'énergie peuvent être identifiées par l'origine de l'arc générateur.

Par dessus le marché, la méthode se prête facilement à la prédiction des spectres à deux quanta. Chaque parallélogramme sur le graphique engendre une transition à deux quanta dont la fréquence correspond à la mi-hauteur du parallélogramme (on doit tourner la figure pour que la direction énergie soit verticale). Le cas A_2X_3 donne naissance à un quadruplet de transitions à deux quanta, avec les intensités relatives 1:4:4:1, en accord avec l'expérience de KAPLAN et MEIBOOM² sur l'alcool éthylique.

Cette méthode s'étend sans difficultés à d'autres systèmes à couplage faible, A_nX_m . On laisse au lecteur le soin de construire (en deux dimensions) les graphiques appropriés aux systèmes $A_nM_pX_q$ (trois groupes de spins inégaux).


Sincères salutations,

Sue

(Mlle) Suzanne Perchery

1. Voir par exemple P. DIEHL, 6 ième Congrès de la R.M.N. Expérimentale (6th E.N.C.) Pittsburgh, 1965; Helv.Chim.Acta, 48, 567 (1965).
2. J. I. KAPLAN et S. MEIBOOM, Phys.Rev. 106, 499 (1957).

Titre: Niveaux d'énergie A_nX_m .



Life Science Systems Process Control Analytical Instruments



System Solutions and Precision Instruments from Bruker:

Bruker NMR/EPR

NMR
MRI
Process Control NMR
EPR

Bruker AXS

Single Crystal Diffraction (SCD)
X-Ray Diffraction (XRD)
X-Ray Fluorescence (XRF)

Bruker Optics

FT-IR
FT-Raman
FT-NIR Quality Control
FT-NIR Life Sciences Systems

Bruker Daltonics

MALDI-TOF & ESI-TOF
FTMS
Ion Trap LC/MS/MS
Proteomics & Genomics Systems

*Innovation for customers
delivered with Integrity*

www.bruker.com

Bruker Worldwide

Bruker is a global group of companies designing, manufacturing and distributing systems for the Life Science, Process Control and Analytical Chemistry markets with production sites in Europe and North America. A network of applications, sales and service offices around the world assures extensive high level customer support.

Our People

are our strength. Through them, Bruker continues to grow and prosper as a leading force in the scientific community. Our scientists, engineers, production staff, and sales team maintain close communication with customers, learning the latest requirements for cutting edge experiments, and incorporating that information into our next generation products. We actively support the exchange of information between our staff and the scientific community through participation in conferences, seminars, symposia, and exhibitions.

Our Solutions

are thorough and complete, concentrating on meeting your analytical requirements. The versatility of Bruker systems and accessories embrace a full range of applications. We offer sophisticated automation products adding high throughput capabilities for industrial applications as well as specialty instruments for unique experiments. Every year brings greater innovations, enhancing the productivity of the labs we serve. Software tools are highly integrated into every product allowing true ease-of-use for both routine as well as research environments.

Our Support

is extensive. The Bruker organization has a well-earned reputation of responsiveness to customers' demands for excellence in every area. Our worldwide organization features offices for sales, application and engineering support in every major area of the world. Whether you need information on new products and accessories, enhancements of existing products, or routine service and support, your local Bruker office is just a quick phone call away.

***Innovation For Customers
Delivered With Integrity***



BRUKER INSTRUMENTS, INC.

44 Manning Road
Billerica, MA. 01821-3991
Phone: (1) 978 667-9580
Fax: (1) 978 667-0985
<http://www.bruker.com>
<http://www.bruker.de>



MASARYK UNIVERSITY
FACULTY OF SCIENCE
 Kotlářská 2, 611 37 Brno
 Czech Republic

Professor Vladimír Sklenář
 Laboratory of Biomolecular Structure
 and Dynamics.
 tel. **420-5-41129579, 41129384
 fax. **420-5-41129506
 e-mail: sklenar@chemi.muni.cz



Prof. B. L. Shapiro
 The NMR Newsletter
 966 Elsinore Court
 Palo Alto, CA 94303
 USA

(received 5/17/99)
 May 17, 1999

HCN TROSY for Nucleic Acids

Dear Barry:

In our previous letter we reported the sensitivity gains obtained with the TROSY¹ approach for ¹H-¹³C correlation in larger oligonucleotides at 500 MHz. The average sensitivity enhancement (TROSY/HSQC) at various sites in a fully ¹³C- and ¹⁵N-labeled ATP binding aptamer containing 40 nucleotides varied between 2.3 and 2.9 at 50°C in the case of aromatic carbons while for C1' sugar carbons the enhancement was substantially smaller (1.3 on average). We did not find the results significant enough to warrant publication, especially since Brutscher and coworkers presented a comprehensive discussion of ¹H-¹³C TROSY correlation for labeled RNA².

Recently, we have decided to resuscitate our TROSY attempts to see the possible merits in applications to ¹³C labeled RNA and DNA oligonucleotides. In continuation of the effort to improve sensitivity of nD HCN and HCNCH experiments for intra-residual correlations of sugar H1' and C1' with N1/9, C6/8 and H6/8 in pyrimidine and purine bases, we have designed their TROSY versions. Previously, the conversion of single- to multiple-quantum evolution in HCN and HCNCH experiments, together with a careful control of homo- and heteronuclear scalar interactions, resulted in up to 400 % increase in signal-to-noise ratio³. The first TROSY results obtained on our spectrometer at 500 MHz were a little bit disappointing. For ATP binding aptamer only pyrimidines showed substantial gains in sensitivity (enhancement factors 2.2 – 2.3) when *out-and-back* TROSY and MQ HCN experiments were evaluated. In the case of purines, the performance of both versions was comparable. For sugar carbons, the performance of the TROSY version of the HCN experiment is always inferior to MQ due to a small value of the chemical shift anisotropy (CSA) of C1'. The sensitivity of through-bond correlating HCNCH experiments, where we combined the MQ evolution on sugars with a TROSY step on bases, was also lower than expected.

Since 11.7 T is far from being optimal for applications of transverse optimized spectroscopy we were wondering what could be achieved at higher fields. Few weeks ago, by courtesy of Bruker Analytik and Wolfgang Bermel, we had a chance to test the experiments at 800 MHz on a AVANCE demo system in Karlsruhe, Germany. As predicted by theoretical calculations (Fig. 1), TROSY approach in both *out-and-back* HCN and *all-the-way-through* HCNCH schemes showed considerable improvement in sensitivity and resolution. Average gain in signal-to-noise ratio was 300 - 400 % for uracils and cytosines and 200-300 % for adenines and guanines when TROSY and MQ versions were compared. The chemical shift labeling in the HCNCH experiment can be arranged as H6/8-to-H1' or H1'-to-H6/8 correlation. Depending on the choice, the TROSY step is placed in the first or second CN evolution interval. The transfer H1'-C1'-N1/9 is designed as the MQ evolution. The S/N gains are not as dramatic as in the case of HCN experiments but average gains of between 22 and 81% were obtained for purine and pyrimidine signals, respectively.

The enhancement factors are not uniform and vary across the RNA structure. To get better understanding of the relaxation properties of C1'(ribose), C6(U and C) and C8(G and A), we have calculated the values of principal components and orientation of CSA using DFT *ab initio* method. The differences in relaxation rates based on the obtained CSA values are in good agreement with the sensitivity relations observed in experimental data. Significant variations observed for individual bases reflect differences of relaxation rates due to local variation of the relaxation environment and of the motional behavior.

Figure 1. Relaxation rates of single quantum coherence (dotted line), multiple-quantum coherence (dashed) and TROSY (solid line) of C6-H6 pair in cytosine (A) and C8-H8 pair in guanine (B) as a function of magnetic field. The calculation included the effect of neighboring protons.




Figure 1. Relaxation rates of single quantum coherence (dotted line), multiple-quantum coherence (dashed) and TROSY (solid line) of C6-H6 pair in cytosine (A) and C8-H8 pair in guanine (B) as a function of magnetic field. The calculation included the effect of neighboring protons.

References

1. K. Pervushin, R. Riek, G. Wieder, and K. Wüthrich, *Proc. Acad. Sci. USA* **94**, 12366 (1997); K. Pervushin, R. Riek, G. Wieder, and K. Wüthrich, *J. Biomol. NMR* **12**, 345 (1998).
2. B. Brutscher, J. Boisbouvier, A. Pardi, D. Marion, and J.-P. Simorre, *J. Amer. Chem. Soc.* **120**, 11845 (1998).
3. R. Fiala, F. Jiang, and V. Sklenář, *J. Biomol. NMR*, **12**, 373-383(1998).

Sincerely,

Radovan Fiala

Rad. Fiala

Jiří Czernek

Jiří Czernek

Vladimír Sklenář

Vladimír Sklenář



BRUKER S.A. - 34, rue de l'Industrie - 67166 WISSENBURG cedex - FRANCE
 Tél. +(33) 03 88 73 68 00 - Télécopieur +(33) 03 88 73 68 79 - Télex 870 639 BRUSPIN
 http://www.bruker.com - e-mail: sadis@bruker.fr

AGENCE RÉGION PARISIENNE : 3, avenue du Général de Gaulle - 91090 LISSES
 Tél. +(33) 01 69 89 85 20 - Télécopieur +(33) 01 60 86 34 90 - e-mail: sarlisse@bruker.fr

WATER INVERSION USING A RADIATION DAMPING CONTROL UNIT

Dear Barry,

(received 4/27/99)

One of the problems frequently encountered when working in aqueous solutions is the efficient inversion of the water signal. This is of particular importance when studying exchangeable protons or bound water molecules in protein samples. The phenomenon of radiation damping will oppose the effect of classical selective excitation pulses (Gaussian, E-BURP...) and will prevent a good excitation of the water signal. The methods that have been proposed in the literature to overcome this problem usually employ techniques that cancel the radiation damping field.

There is, however, another way to achieve a clean water inversion which is based on the excitation of the water signal by its own radiation damping field. The idea, originally proposed by Abergel et al.¹, uses an electronic circuit to detect and manipulate the radiation damping field. The net result is that the resulting field will not bring the water magnetization along the +z axis (classical radiation damping case) but along the -z axis (water inversion case).

The results we have obtained using our newly developed radiation damping control unit are quite impressive and lead to a very good water inversion. When taking the difference between two spectra recorded with the water magnetization along the +z or the -z axis, only the protons in interaction (chemical exchange or dipolar interaction) with water molecules appear in the resulting spectrum. The spectra we have obtained on a 10 mM Cytochrome c sample² (Figure 1) exhibit a signal to noise superior to the one obtained using alternative techniques^{3,4}.

O. BORNET, F. GUERLESQUIN, C. BREVARD, M. PIOTTO

¹ D. Abergel, A. Louis-Joseph and J-Y Lallemand, *Chem. Phys. Lett.*, **262**, 465 (1996).

² O. Bornet, F. Guerlesquin and M. Piotto, *J. Magn. Reson.*, in press.

³ C. Dalvit and U. Hommel, *J. Magn. Reson. B*, **109**, 334 (1995).

⁴ C. Dalvit, *J. Magn. Reson. B*, **112**, 282 (1996).






Figure 1: Comparison of the results obtained on a 10 mM sample of horse heart ferrocyanochrome *c* in 90% H₂O and 10% D₂O at 300 K. A 500 MHz Avance Bruker spectrometer equipped with a HCN probe and self shielded triple axis gradients was used. A) Reference 1D spectrum recorded with 8 scans and presaturation. B) 1D NOE spectrum recorded using a selective spin echo for water inversion^{3,4} with 64 scans and a 100 ms NOE mixing time. A 50 ms Gaussian pulse was used to select the water magnetization. C) 1D NOE spectrum recorded using a radiation damping control unit for water inversion with 64 scans and a 100 ms NOE mixing time.




Intuitive interface for all kinds of processing, searching & databasing



Tile feature in Spectrum Window lets you view several spectra at once

Standardize data input with Data Forms Manager



Comprehensive Database Window

ACD/NMR Manager

Take raw NMR data directly from the spectrometer and process it using a wide array of tools!

- Import the major FID and FT formats, including Varian, Bruker, JEOL, Chemagnetics, WinNMR, Nuts, Lybrics, JCAMP, etc.
- Use zero filling, weighting functions, Fourier Transform.
- Apply automated or manual phase and baseline correction, peak referencing.
- Automated or manual peak picking and integration.
- Manual or automatic peak assignment.
- Annotate the spectrum using either peak or region selection, offering unique capabilities to assign and annotate polymer spectra and other complex spectral curves.
- Create macro programs for completely automated data processing.
- Dual Window display of calculated and experimental spectra.
- **Plus** - All new capabilities of NMR Processor

Innovative storage and management software for all kinds of experimental data!

Create a *comprehensive* database of experimental spectra:

- Store different kinds of spectra in the same data base;
- **New!** - Assign several spectra to a single data base entry;
- Establish over 16,000 fully searchable user-defined fields;
- Protect databases with view-only and read-write passwords;
- Search database according to:
 - spectrum, sub-spectrum, spectral parameters;
 - molecular structure or sub-structure;
 - formula, molecular weight and user-defined fields; and any special fields corresponding to spectrum type
- Search through several databases at a time.
- Generate and save lists from searches, perform various manipulations with lists (duplication, merging, subtraction, intersection) to make complex searches.
- **New!** - Standardize and simplify record entry with the Data Forms Manager.

MS, UV-IR, 2D-NMR Modules also available.




NMR

IR

UV-Vis

MS

Raman





Visionary Software for Scientists

Advanced Chemistry Development


NMR

ACD/NMR Processor

Import almost any 1D FID format




Multiple window display includes Synchronize Axes button




Reference peaks table

Name	Offset	Multiplicity
CHCl3-d4	7.25	1
Acetic acid-d4	2.33	3
Acetone-d6	2.14	3
Acetonitrile-d3	1.94	3
Dimethyl sulfoxide-d6	2.50	3
Methanol-d4	3.30	3


Intuitive interface for processing



Rectangular zooming



Create professional reports



Take raw NMR data directly from the spectrometer and process it using a wide array of tools!

- Import the major FID and FT formats: Varian, Bruker, JEOL (EX, GX, Alpha, Lambda), Felix, Chemagnetics, WinNMR, Nuts, Lybrics, JCAMP, GE Omega, Tecmag, and ASCII.
- Apply automated or manual phase and baseline correction.
- Automated or manual peak picking, referencing, & integration.
- Fast click-and-drag assignment between structure and spectrum.
- Annotate the spectrum using either peak or region selection, offering unique capabilities to assign and annotate polymer spectra and other complex spectral curves.
- Create macro programs for completely automated data processing.
- **New!** - Interactive apodization takes the guesswork out of processing.
- **New!** - Five varieties of weighting functions.
- **New!** - Multiple window display.
- **New!** - Add text in Comments field.
- **New!** - Synchronized zoom over multiple windows.

SpecViewer
www.acdlabs.com
FREE!

Advanced Chemistry Development
T: 416 368-3435
F: 416 368-5596
Toll Free: 1 800 304-3988
info@acdlabs.com
www.acdlabs.com

Raman



Department of Chemistry
Professor Gideon Fraenkel
 Phone: 614-292-4210
 e-mail: fraenkel@mps.ohio-state.edu

April 28, 1999
 (received 5/5/99)


Newman and Wolfrom Laboratory
 100 West 18th Avenue
 Columbus, OH 43210-1185
 Phone 614-292-2251
 FAX 614-292-1685
 TELEX 332911
 Answer Back Code: OSU CHEM UD

Dr. B. L. Shapiro
 The NMR Newsletter
 966 Eismore Court
 Palo Alto, CA 94303-3410


Dear Barry:

Solvated allylic lithium compounds do not ordinarily exhibit one bond ^{13}C , ^6Li (or ^7Li) coupling. With coordinated lithium sited normal to the allyl plane, the "s" contribution to any C, Li covalency would be minimized, hence spin coupling is too small to be detected. The exceptions which show ^{13}C lithium coupling are the few allylic lithium compounds within which lithium is internally coordinated to a pendant ligand, see **1**. In these latter cases we proposed¹ (confirmed with X-ray crystallography²) that restricted stereochemistry of lithium coordination favors lithium to be somewhat off the vertical axis (40°) normal to the allyl plane. This increases the "s" character associated with a small C, Li covalence, allowing for detectable (3 Hz) ^{13}C , ^6Li coupling.

By lengthening the ligand tether by one carbon, we now have a system in which there is ^{13}C , ^7Li coupling to both termini of an allylic lithium compound. It is produced by metalating **2** with CH_3Li at CH_2Si .




1



2

The lithium compound is a monomer in solution (freezing point) and in the solid (X-ray). The latter shows a covalent species with Li bonded to the terminus with silyl endo. In solution, ^7Li is coupled to the two ^{13}C termini differently, 4.3 Hz and 4.9 Hz, up to 300 K. The terminal ^{13}C shifts of 64.88 and 68.38 don't look like a localized system. However, all can be explained.

In contrast to the solid state, in solution there is an equilibrium mixture of two covalent species in similar concentration rapidly interconverting via a lithium 1,3-sigmatropic shift, see $3a \rightleftharpoons 3a'$, $3b \rightleftharpoons 3b'$. Relative to the NMR time scale bimolecular C, Li exchange is slow as is also rotation around the $\text{C}_1\text{-C}_2$ and $\text{C}_2\text{-C}_3$ bonds. Our sample has ^{13}C in natural abundance. Thus, the two pairs of separately interconverting isotopomers are $3a \rightleftharpoons 3a'$, $3b \rightleftharpoons 3b'$. So each observed coupling constant is the average over two interconverting species, $3a$, a' and $3b$, b' . They are different because the bond rotation rates which convert $3a$ to $3b$ or $3a'$ to $3b'$ are slow in this sample. The observed shifts are likewise averages over $3a$, a' and $3b$, b' .



All of this comes from variations in the stereochemistry of lithium solvation.

Best regards,

Gideon

Gideon Fraenkel
Professor of Chemistry

Roland

Roland Fleischer
Research Associate

1. Fraenkel, G., Qiu, F. *J. Am. Chem. Soc.* 1997, 119 3571-3579.
2. Watch your JACS this summer.

Thomas L. James, Ph.D.
Chairman
Department of Pharmaceutical Chemistry
University of California, San Francisco
926 Medical Sciences Building
San Francisco, CA 94143-0446
(415) 476-1916
FAX: (415) 502-4690
ELECTRONIC MAIL: james@picasso.ucsf.edu
Webpage: <http://picasso.ucsf.edu:80/james/>

UCSF
University of California, San Francisco ... A Health Science Campus

May 27, 1999

Dr. B.L. Shapiro
The NMR Newsletter
966 Elsinore Court
Palo Alto, CA 94303

Re: For Sale: 4.7 Tesla, 330 mm Horizontal Clear-Bore Magnet & Accessories

Dear Barry:

I am nearing the final step in my multiyear ease out of the *in vivo* NMR business. As a consequence, I have for sale a horizontal bore instrument which is adequate for studies with animals up to the size of young rabbits. While the console (Nalorac Quest) is aged and probably of little value to others, the Oxford magnet is really excellent. The magnet has actively shielded shim coils and gradient coils (15 G/cm maximum), and it has a nice "sweet" spot according to users. The magnet is currently operating, but we will probably decommission it within a month. Other components, such as frequency synthesizer and power amplifiers, may be useful to some people as well. More details about the available equipment and contacts regarding procurement can be found at the web site <http://picasso.ucsf.edu/magnetforsale.pdf>.

Sincerely yours,

Tom

Thomas L. James, Ph.D.
Chairman, Department of Pharmaceutical Chemistry
Professor of Chemistry, Pharmaceutical
Chemistry and Radiology

SIMON FRASER UNIVERSITY

INSTITUTE OF MOLECULAR
BIOLOGY AND BIOCHEMISTRY



8888 UNIVERSITY DRIVE
BURNABY, BRITISH COLUMBIA
CANADA V5A 1S6
Telephone: (604) 291-5630
Fax: (604) 291-5583

April 12, 1999
(received 4/27/99)

Dr. B. L. Shapiro
The NMR Newsletter
966 Elsinore Court

Spectral Estimation of NMR Relaxation

Dear Dr. Shapiro

If a monoexponential decay has the following form:

$$f(t) = Ae^{-\lambda t} + B$$

then the one sided continuous Fourier Integral Transform gives the transform as:

$$F(\omega) = A(\lambda - i\omega) / (\lambda^2 + \omega^2) + 2\pi B\delta(\omega)$$

and hence the decay rate λ can be extracted as:

$$\lambda = -\omega \Re(F(\omega)) / \Im(F(\omega))$$

for any nonzero angular frequency ω .

In "Exponential Analysis of physical phenomenon" by Andrei A. Istratov and Oleg F. Vyvenko, in Review of Scientific Instruments, Vol. 70, No. 2, Feb. 1999, p1233-1257, (with 344 references), the authors advocate this method as the best method for determination of a monoexponential decay rate. This relationship is a very good approximation where the discrete Fourier Transform is used with long records of 256 points but is not very approximate when short records like four, eight, twelve points are collected, like in NMR relaxation.

In our opinion, current wisdom in NMR is represented by the ideas presented in "Optimization of Magnetization Transfer Experiments to Measure First-Order Rate Constants and Spin-Lattice Relaxation Times", by H. Katki et al. in NMR in Biomedicine, 9, 135-139 (1996), and "The Effects of Imperfect Saturation in Saturation Recovery T1 Measurements" by A. Roscher et al. in J. Mag. Res. 118A p108-112 (1996).

However, it can be shown that for a short record DFT, the decay rate can be extracted from data as:

$$\lambda = \ln(-\Re(F(\omega)) / \Im(F(\omega)))$$

where, ω is selected to be at the half Nyquist angular frequency. The DFT can be thought of as a matrix of complex values and hence this latter expression can be used to construct a more rational expression for the extraction of λ from data. If the decay transient is measured at points with equal time increments and if the time increment is taken as the unit of time, λ can be extracted from data as follows:


For four points, f_1, \dots, f_4

for twelve point, f_1, \dots, f_{12}

$$\lambda = -\ln((f_2 - f_4 + f_6 - f_8 + f_{10} - f_{12}) / (f_1 - f_3 + f_5 - f_7 + f_9 - f_{11}))$$

These formulae provide a non iterative estimate of λ which is independent of any assumption concerning the error distribution within the data.

Sincerely,



David Naugler
dnaugler@sfu.ca



Robert J. Cushley
cushley@sfu.ca

Reply by:

Dr. Andrei A. Istratov,
Institute of Physics of St. Petersburg State University
Ulianovskaya 1
St. Petersburg, 198904
Russia

Re: Spectral Estimation of NMR Relaxation

David,

This is indeed a good point. Equation of the form $\lambda = -\ln((f_2 - f_4) / (f_1 - f_3))$ appears in Sec. VI.A.2 of the review, and the obvious disadvantage of this formula is that it uses only four points, no matter how many are actually measured. Equation of the type $\lambda = -\ln((f_2 - f_4 + f_6 - f_8 + f_{10} - f_{12}) / (f_1 - f_3 + f_5 - f_7 + f_9 - f_{11}))$ is much better in this respect, and I agree that it deserved to be included in the review. Thank you for pointing this out. You know, after I read your e-mail I realized that we should have discussed in more detail in the review what to do if the number of experimental data points is limited to, e.g., 6, 8, or 10. This is a very special problem, and only few algorithms can be used.

Another point that was discussed in Sec. VI (page 1242, left column, at the top of the page), but was not emphasized in the discussion, is that all methods for multiexponential analysis are suitable for monoexponential analysis. For this reason, such a common technique as exponential nonlinear least squares fit would be a good algorithm to analyze a set of 8 or 10 data points. But it is surely not as fast as your formula.

Thank you for your comments,

Sincerely,



Andrei Istratov
istratov@socrates.berkeley.edu

New Stable Isotope Products

Chemists in the Stable Isotopes group at Aldrich are committed to serving the product needs of our customers throughout the spectrum of scientific fields. We are continually expanding our product line to satisfy all of your research needs, whether they be in deuterated solvents for NMR spectroscopy or labeled specialty chemicals. Below is a sampling of the most recent additions to this product line. If you would like further information or have a need for a labeled material that is not currently available, please call us at **800-231-8327 (USA)** or visit our Web site at **www.sigma-aldrich.com**.

Cat. No.	Isotopic Purity Atom % D	NMR Reference Standard	Tube Dimension	Pkg Size	\$
48,775-9*	99.9	5% chloroform in acetone- d_6	5mm x 8in.	1 ea (0.70mL)	66.90
48,717-1*	99.9	1% chloroform in acetone- d_6	8mm x 8in.	1 ea (1.95mL)	128.70
48,716-3*	99.9	1% chloroform in acetone- d_6	5mm x 8in.	1 ea (0.70mL)	66.90
48,714-7*	99.9	1% 1,2-dichlorobenzene in acetone- d_6	5mm x 8in.	1 ea (0.70mL)	77.20
48,712-0*	99.8	0.1% ethylbenzene and 0.01% tetramethylsilane in chloroform- d	8mm x 8in.	1 ea (1.95mL)	149.30
48,710-4*	99.8	0.1% ethylbenzene and 0.01% tetramethylsilane in chloroform- d	5mm x 8in.	1 ea (0.70mL)	92.70
48,711-2*	99.8	0.1% ethylbenzene and 0.01% tetramethylsilane in chloroform- d	3mm x 8in.	1 ea (0.23mL)	108.10
48,713-9*	99.9	0.1 mg/mL gadolinium(III) chloride, 0.1% DSS and 1% water in deuterium oxide	5mm x 8in.	1 ea (0.70mL)	154.50
48,715-5*	99	1% iodomethane- ^{13}C , 1% trimethyl phosphite and 0.2% chromium(III) acetylacetonate in chloroform- d	3mm x 8in.	1 ea (0.23mL)	145.00

Cat. No.	Isotopic Purity Atom % D	Name	Pkg Size	Pkg	\$
32,992-4		Tetrakis(trimethylsilyl)silane, 98% (solid-state NMR reference standard)	1g 10g	screw-cap bottle screw-cap bottle	18.30 104.10
32,673-9		Hexamethyldisiloxane, 99.5+%, NMR grade	25g 100g	screw-cap bottle screw-cap bottle	12.00 33.00
T2,400-7		Tetramethylsilane, 99.9+%, NMR grade, A.C.S. reagent	25g 100g	screw-cap bottle screw-cap bottle	19.90 56.30
17,883-7		3-(Trimethylsilyl)-1-propanesulfonic acid, sodium salt, 97%	1g 5g	screw-cap bottle screw-cap bottle	17.30 67.10
18,033-5		3-(Trimethylsilyl)propionic acid, sodium salt, 99+%	1g 5g	screw-cap bottle screw-cap bottle	34.10 126.90
26,991-3	98	3-(Trimethylsilyl)propionic-2,2,3,3- d_4 acid, sodium salt (TSP- d_4)	1g 5 x 1g	screw-cap bottle screw-cap bottle	37.85 120.20

* Manufactured by **ISOTECH**.



NEW Aldrich Deuterated NMR Solvents Brochure

- ◆ Complete listing of all deuterated solvents, including those with 100.0 atom % D
- ◆ Broadened product line of deuterated products for NMR, with many new products and package sizes listed
- ◆ Reference spectra and standards

To receive your copy of the **New Deuterated Solvents for NMR** brochure (available in July 1999), please e-mail dszydel@sial.com or call 800-771-6737 (USA) ext. 5784.



ALDRICH®

chemists helping chemists in research & industry

P.O. Box 355, Milwaukee, WI 53201 USA Telephone: 414-273-3850 • 800-558-9160 Fax: 414-273-4979 • 800-962-9591 Web Site: www.sigma-aldrich.com

Aldrich is a member of the Sigma-Aldrich family.

Sure/Seal™ bottles for NMR solvents

Aldrich now has our most popular NMR solvents available in Sure/Seal™ bottles. Our well-known, exclusive packaging system safeguards against leakage and exposure to moisture and oxygen. In addition, it virtually eliminates the loss of TMS in our chloroform listings. The Sure/Seal™ bottle system allows for easy storage, handling, and dispensing of product using standard syringe transfer techniques. You are assured of a quality product, from the first time you open the package to the final product usage.

To cover all your needs, we offer all of our NMR solvents in a number of other packaging options as well, including screw-cap bottles and individual ampoules. For more information and a complete list of our deuterated solvents, please see our 1998-99 Catalog/Handbook of Fine Chemicals, call us at 800-231-8327 (USA), or visit us on the Web at www.sigma-aldrich.com.



Please call us today to get a competitive quotation for your annual NMR solvents requirements.

Cat. No.	Name, Isotopic Purity	Unit Size and Price†
44,486-3	Acetone- d_6 , 99.9 atom % D	100mL \$175.10
15,179-3	Acetone- d_6 , 99.5 atom % D	100g \$146.20
15,180-7	Acetonitrile- d_3 , 99.6 atom % D	100g \$228.60
15,181-5	Benzene- d_6 , 99.6 atom % D	100g \$175.10
15,182-3	Chloroform- d , 99.8 atom % D	150g \$31.90
22,578-9	Chloroform- d , 99.8 atom % D (contains 0.03% v/v TMS)	150g \$31.90
43,487-6	Chloroform- d , 99.8 atom % D (contains 0.1% v/v TMS)	150g \$31.90
15,183-1	Chloroform- d , 99.8 atom % D (contains 1.0% v/v TMS)	150g \$31.90
15,189-0	Deuterium oxide, 100.0 atom % D	125g \$169.90
15,188-2	Deuterium oxide, 99.9 atom % D	125g \$72.10
15,194-7	Methyl- d_3 alcohol- d , 99.8 atom % D	100g \$451.10
15,187-4	(Methyl sulfoxide)- d_6 , 99.9 atom % D	100g \$142.10
17,594-3	(Methyl sulfoxide)- d_6 , 99.5+ atom % D	100g \$133.90

† This is the quantity packaged under nitrogen in the Sure/Seal™ bottle. We offer many other units and packaging options. Please contact us for more details.

Sure/Seal is a trademark of Sigma-Aldrich Co.

Teflon is a registered trademark of E.I. du Pont de Nemours & Co., Inc.



ALDRICH®

chemists helping chemists in research & industry

P.O. Box 355, Milwaukee, WI 53201 USA Telephone: 414-273-3850 • 800-558-9160 Fax: 414-273-4979 • 800-962-9591 Web Site: www.sigma-aldrich.com

Aldrich is a member of the Sigma-Aldrich family.



UNIVERSITY OF ALBERTA

Dr. B. L. Shapiro
The NMR Newsletter
966 Elsinore Court
Palo Alto, CA 94303
USA

May 14, 1999
(received 5/20/99)

Temperature Calibration of a Solid State NMR probe

Dear Barry,

We have recently started using ^2H solid state NMR at low temperatures and wanted to calibrate the temperature in the solids probe. Typically, for high resolution work, one uses methanol and measures the difference between the OH and CH_3 resonances. Since our solids probe is ^2H , we were interested in how accurate this procedure would be using the OD and CD_3 resonances.

The central issue in transferring ^1H temperature calibration curves to ^2H NMR relates to whether $\gamma_{\text{H}} / \gamma_{\text{D}}$ remains constant as a function of temperature. If it does remain constant, our calibration of the probe temperature is simple. We expect a value for $\gamma_{\text{H}} / \gamma_{\text{D}}$ of approximately 6.51. Having carefully calibrated our ^1H -indirect detection probe using methanol, we ran a spectrum of Val- H_8 at different temperatures and measured the difference between the HDO peak and the center of the Val- $(\text{CH}_3)_2$ resonances. Our calculations indicated that $\Delta\delta = -3.73T + 1223.8$ for the indirect detection probe (Figure 1A). Next, we measured the difference between the HDO peak and Val- $(\text{CD}_3)_2$ resonances in the solids probe. We found $\Delta\delta = -0.55T + 194.0$ (Figure 1B). The slope and y-intercept of the ^2H curve is different from the ^1H curve by an average factor of 6.5. This information indicates to us that $\gamma_{\text{H}} / \gamma_{\text{D}}$ is not a function of temperature. Therefore, we can use the calibration of the indirect detection probe to get an accurate temperature inside our solids probe.

Department of Biochemistry

474 Medical Sciences Building • University of Alberta • Edmonton • Canada • T6G 2H7

Telephone: (780)

• Fax: (780) 492-0886

www.ualberta.ca




Figure 1: Temperature calibration curves for (A) ^1H indirect detection probe and (B) solid state ^2H probe. Shown is the difference in chemical shift between the HDO peak and the center of either Val-(CH_3) $_2$ (^1H indirect detection probe) or Val-(CD_3) $_2$ (^2H probe) versus temperature. The calibrated temperature inside the ^1H indirect detection probe is shown in (A) whereas the temperature reading on the VT console was used on the solids probe in (B).

Sincerely,

Carolyn Slupsky

Brian Sykes

Water Research Institute

Sengen 2-1-6
Tsukuba
Ibaraki 305-0047
Japan
Ph: (81-298) 58 6186
FAX: (81-298) 58 6144
e-mail: wprice@wri.co.jp
May 12, 1999
(received 5/19/99)

Dr. B. L. Shapiro
The NMR Newsletter
966 Elsinore Court
Palo Alto, CA 94303
USA

PGSE NMR Gradient Pulses and Artifactual Diffraction Peaks

Dear Barry,

When doing diffusion measurements especially with very high amplitude gradient pulses, phase disturbances are often noted. One commonly used means of 'solving' this problem is to hide it by computing absolute value spectra. It is interesting to consider this situation and this 'solution' a little further. For example, if we perform an analysis based on the short gradient pulse approximation and include the effects of a phase-shift, ϕ , due to the effects of a gradient mismatch, Δq , and sample movement, Δr , between the first and second gradient pulses in the PGSE sequence we get (1)

$$E(\mathbf{q}, \Delta) = \int \rho(\mathbf{r}_0) \int P(\mathbf{R}, \Delta) e^{i[2\pi \mathbf{q} \cdot \mathbf{R} + \phi]} d\mathbf{R} d\mathbf{r}_0 \quad [1]$$

where $\rho(\mathbf{r}_0)$ is the spin density, \mathbf{R} is the dynamic displacement defined by $\mathbf{r}_1 - \mathbf{r}_0$ (the starting and finishing positions of a spin with respect to the first and second gradient pulses), $P(\mathbf{R}, \Delta)$ is the average propagator (i.e., the probability that a spin will move by a displacement \mathbf{R} during the time interval Δ) and the phase term can be expressed as

$$2\pi \mathbf{q} \cdot \mathbf{R} + \phi = 2\pi [(\mathbf{q} + \Delta \mathbf{q}) \cdot (\mathbf{r}_0 + \mathbf{R} + \Delta \mathbf{r}) - \mathbf{q} \cdot \mathbf{r}_0]. \quad [2]$$

Taking the gradient to be oriented along the z -direction and assuming that Δq is a magnitude mismatch (i.e., parallel to \mathbf{q}), Eq. [1] can be re-written as

$$E(q, \Delta) = \underbrace{\int P(Z, \Delta) e^{i2\pi q Z} dZ}_{E_{\text{diff}}(q, \Delta)} \underbrace{\left\{ e^{i2\pi(q + \Delta q)\Delta z} \right\}}_{\phi(\Delta z)} \underbrace{\int \rho(z_1) e^{i2\pi \Delta q z_1} dz_1}_{E_{\text{phase}}(\Delta q)}. \quad [3]$$

The first term is the attenuation due to diffusion, the second term is the residual (net) phase-shift due to sample movement (which can be removed by computing the absolute value spectrum), and the third term is the residual phase-twist resulting from the gradient pulse mismatch. For a cylindrical sample of length l centred in the gradient we obtain

$$E_{\text{phase}}(\Delta q) = \text{sinc}(\pi \Delta q l). \quad [4]$$

Some example plots of $\ln(|E(q, \Delta)|)$ as they would appear if they were determined from absolute value spectra are shown in Figure 1. It has been simplistically assumed that Δq is a fixed proportion of q . It can be clearly seen that as the mismatch increases, the measured (i.e., apparent) diffusion coefficient increases (i.e., compare the slopes at small q values).




Figure 1 Plots of $\ln(|E(q, \Delta)|)$ calculated using Eqs. [3] and [4] with $D = 1 \times 10^{-14} \text{ m}^2\text{s}^{-1}$, $\Delta = 20 \text{ ms}$, $l = 5 \text{ mm}$, $\Delta z = 0 \text{ mm}$ and $\Delta q = 0$ (—) (i.e., no mismatch, this corresponds to $E_{\text{diff}}(q, \Delta)$), $\Delta q = q \times 2 \times 10^{-4}$ (····) and $\Delta q = q \times 1 \times 10^{-3}$ (----). With such a low value of D the diffusive attenuation should be almost negligible as shown for the case of $\Delta q = 0$, however, the residual phase-shift resulting from gradient pulse mismatch causes enormous artifactual attenuation resulting in a huge increase in the apparent diffusion coefficient (i.e., if the attenuation data is interpreted assuming that all of the observed attenuation resulted from $E_{\text{diff}}(q, \Delta)$). At larger q values (and therefore Δq) artifactual diffraction peaks are evident. The degree of mismatch strongly influences the position of the diffraction peaks with the diffractive peaks appearing at smaller values of q as the degree of mismatch increases.

We note that the $E_{\text{phase}}(\Delta q)$ term is not removed by computing absolute value spectra. However, if absolute value spectra are used the negative excursions in the attenuation plot become positive. This results in an attenuation plot that contains sharp minima similar to the diffractive minima that can be observed for species undergoing restricted diffusion (2). Thus, it is important to observe the PGSE spectra in phase-sensitive mode to distinguish between real and artifactual diffraction peaks in PGSE attenuation plots.

References

1. P. T. Callaghan, PGSE_MASSEY, a sequence for overcoming phase instability in very-high-gradient spin-echo NMR, *J. Magn. Reson.* **88**, 493-500 (1990).
2. W. S. Price, K. Hayamizu, H. Ide, and Y. Arata, Strategies for diagnosing and alleviating artifactual attenuation associated with large gradient pulses in PGSE NMR diffusion measurements, *J. Magn. Reson.* (1999) In press.

Please credit this to the account of Prof. Y. Arata.

Yours sincerely

W. S. Price

William S. Price

The University of Virginia

ROBERT G. BRYANT
Commonwealth Professor

Department of Chemistry
McCormick Road
Charlottesville, VA 22901
(804) 924-1494
rgb4g@virginia.edu


May 4, 1999
(received 5/10/99)

Barry L. Shapiro
Editor, The NMR Newsletter
966 Elsinore Court
Palo Alto, CA 94303

Oxygen T_{1e} from Magnetic Relaxation Dispersion

Dear Barry,

Oxygen contributions to magnetic relaxation are ubiquitous and very well known. Perhaps not so well known is the magnetic field dependence of the oxygen contribution to relaxation rates. For small molecule systems, the correlation time for the electron-nuclear coupling will become the electron T_{1e} in the case that the electron relaxation rate is rapid compared with correlation times for translational diffusion. Our new MRD spectrometer permits an easy measurement that reports the T_{1e} as the inflection point in a relaxation dispersion curve. Data for oxygen contributions in cyclohexane are shown in the Figure below for air and a sample prepared at 1 atm of oxygen. The oxygen g-factor¹ is 2.0021, so that the top axis of the figure may be labeled easily with the electron Larmor frequency, ω_s . On this axis, the inflection point corresponds to a correlation time of 5 ps, which is much shorter than translational correlation times for either the liquid or the dioxygen molecules in this sample. Thus, this inflection provides a convenient measurement of the oxygen T_{1e} . We have observed similar results in other solvents, such as water, methanol, etc.



Ching-Ling Teng
Ching-Ling Teng

Bob
Robert G. Bryant

Shawn Wagner
Shawn Wagner

1. T. J. Cook, B. R. Zegarski, W. H. Breckenridge, T. A. Miller, J. Chem. Phys. 58, 1548 (1973).



CENTRE NATIONAL
DE LA RECHERCHE
SCIENTIFIQUE



UNIVERSITE CLAUDE-BERNARD LYON 1



LABORATOIRE DE RESONANCE MAGNETIQUE NUCLEAIRE

METHODOLOGIE ET INSTRUMENTATION EN BIOPHYSIQUE - CNRS UPRESA 5012

Villeurbanne le 20 Avril 1999

(received 4/27/99)

Dr B. L. SHAPIRO

The NMR Newsletter

966 Elsinore Court

Palo Alto CA 94 303

Pseudorandom k -space sampling

Dear Dr. Shapiro,

NMR imaging is particularly sensitive to motion effects, especially functional imaging. Random gross motions of the patient may induce pseudorandom sampling when using cartesian scanning. In the present analysis, we assume that random motions do not occur during sampling along a trajectory but only during the time slots between successive k -space trajectories. Figure 1 shows the cells computed with the Voronoi algorithm (1) for a pseudorandom sampling distribution. The cell areas represent the inverse of the sampling density. In this example, random rotations, occurring in the time slots between successive k -space horizontal trajectories were only considered.




Figure 1 : Voronoi cells associated with each sample of a pseudorandom k -space sampling distribution obtained in presence of rotation effects. In the absence of such rotation effects, the samples would be on a Cartesian grid, all the cells would be rectangular and would have the same area.

To reconstruct the correct image, the Bayesian estimator (2) was found more efficient than the gridding algorithm (3). A comparison between these two methods will be discussed during the defense of Marc Bourgeois' thesis on June 9.

Yours Sincerely,

Marc Bourgeois

André Briguet

Danielle Graveron-Demilly

- (1) M. Bourgeois, Y. Crémillieux, C. Delon-Martin et al, International workshop on Sampling Theory and Applications, pages 34-36, Aveiro, Portugal June 1997
- (2) F.T.A.W. Wajer, G.H.L.A. Stijman, R. de Beer et al, in Proc ASCL, pages 147-153, 1998
- (3) J.D.O. Sullivan, IEEE Trans. Med. Imag. 4, 200-207, 1985.



TEXAS A&M UNIVERSITY

Department of Chemistry
College Station, Texas 77843-3255

(409) 845-2011

FAX (409) 845-4719

May 14, 1999 (received 5/21/99)

Glycine Revisited

Dear Dr. Shapiro,

Glycine ($C_2H_5NO_2$) is perhaps the simplest of all amino acids. There is no center of chirality. While it has many uses in other areas, glycine is a useful and convenient sample for solid state NMR. In particular, it can be used to check and calibrate spectrometer performance or, with the various isotopic enrichments commercially available, to quickly implement and demonstrate a variety of different techniques, such as rotational resonance, REDOR, triple resonance experiments, etc.

Glycine is known to exist in three forms (α , β , and γ) in the solid state. These have been studied by both x-ray(1,2,3) and neutron diffraction(4). The α form has hydrogen bonded double layers of molecules which are packed by van der Waals forces. The β form is unstable. It has single molecular layers, whose internal arrangements are the same as in the α form. The γ form crystalizes with a trigonal hemihedral symmetry. The γ form is stable at room temperature but will convert to the α form upon heating to near its decomposition temperature of 165 C. The unstable β form is reported to change to α upon heating, by mechanical shock, or by grinding. The three forms can also be identified by the N-H stretching frequency in the infrared spectra(2).

Figure 1 shows the ^{13}C CP/MAS spectrum of glycine (as is). This sample, of unknown age and handling, originally came from Aldrich Chemical Company. The spectrum shows two carbonyl resonances. The α form is assigned to the resonance at 176.1 ppm while the γ form is assigned to 174.3 ppm. The methylene resonances overlap but may be separated by using the differences in the proton spin lattice relaxation rates. Relaxation time constants are given in Table 1.

Relatively pure samples of either the α or γ form are easy to prepare from such a glycine mixture. The α form can simply be recrystallized by slow evaporation from aqueous solution while the γ form is rapidly recrystallized from aqueous solution with the addition of ethanol. The literature values for the N-H stretching frequency in the infrared spectra have been used to identify the purified forms. The N-H stretch appears at 3173.4 wavenumbers for α glycine and at 3107.8 wavenumbers for γ glycine. The IR spectrum of each form is shown in Figure 2. The above ^{13}C spectral assignments were made based upon the identification of these purified samples.

The NMR relaxation parameters differ significantly between the α and γ forms. However, the simple preparation of either allows the spectroscopist a choice of which to use.

Sincerely,

Robert E. Taylor
Research Instrumentation Specialist

Please credit this contribution to the subscription of Steve Silber.




Figure 1: ^{13}C CP/MAS spectrum of glycine (as is) from Aldrich Chemical Company.




Figure 2: The infrared spectra using KBr pellets of γ glycine (top) and α glycine (bottom).

Table 1: Relaxation Times

Glycine	^1H T1	^1H T1 rho	^{13}C T1	^{13}C T1 rho
α	0.21 sec	45 msec	12 sec (C=O) 4.6 sec (NCH ₂)	220 msec (C=O) 37 msec (NCH ₂)
β	3.5 sec	5.7 msec	43 sec (C=O) 20 sec (NCH ₂)	300 msec (C=O) 37 msec (NCH ₂)

All NMR measurements were made with a Bruker MSL-300 spectrometer with ^1H and ^{13}C ninety degree pulsewidths of 4 usec.

References:

1. Richard E. March, *Acta Crystallogr.* (1958), **11**, 654.
2. Y. Iitaka, *Acta Crystallogr.* (1960), **13**, 35.
3. Y. Iitaka, *Acta Crystallogr.* (1958), **11**, 225.
4. A. Kvik, W. M. Canning, T. F. Koetzle, and G. J. B. Williams, *Acta Crystallogr.* (1980), **B36**, 115.

**Address all Newsletter
correspondence to:**

Dr. B. L. Shapiro
The NMR Newsletter
966 Elsinore Court
Palo Alto, CA 94303.
650-493-5971* - Please call
only between 8:00 am and
10:00 pm, Pacific Coast time.

Deadline Dates

No. 490 (July) 25 June 1999
No. 491 (Aug.) 23 July 1999
No. 492 (Sept.) 20 Aug. 1999
No. 493 (Oct.) 24 Sept. 1999
No. 494 (Nov.) 22 Oct. 1999

* Fax: 650-493-1348, at any hour. Do not use fax for technical contributions to the Newsletter, for the received fax quality is very inadequate.

* E-mail: shapiro@nmrnewsletter.com




The Newsletter's fiscal viability depends very heavily on the funds provided by our Advertisers and Sponsors. Please do whatever you can to let them know that their support is noted and appreciated.

Mailing Label Adornment: Is Your Dot Red ?

If the mailing label on your envelope is adorned with a large **red dot**: this decoration means that you will not be mailed any more issues until a technical contribution has been received.

JEOL Can Give You the Data You Need From Your Desktop PC or MAC



The **Eclipse+** NMR Spectrometer can be operated anywhere there is a computer on the local network. The **Single Window Automation** pictured above can be used with a single mouse click to select the sample from the auto-sample changer, gradient shim on any probe, run the selected experiment, and plot the data on any network postscript printer. Need more data, click another button and the **Eclipse+** is off to do your work - and you have not left your office. Contact us at nmr@jeol.com or visit our web site at www.jeol.com.

JEOL USA, Inc., 11 Dearborn Road, Peabody, MA 01960
Tel: 978-535-5900 Fax: 978-536-2205
email: nmr@jeol.com www.jeol.com

JEOL